

**Hong Kong Institute of Vocational Education (Tsing Yi)
Department of Information Technology**

Higher Diploma in Software Engineering (IT114105)

Final Year Project (ITE4116M)

Final Report

Title: Remote Health Monitoring System

Supervisor: **CHEUNG KWOK FAI
BILLY**

Co-supervisor: Ms Chiu

Student: Wong Ho Shun Sang
 CUI Tianlang
 Yau Yuk ming
 Chan Hing Chak

We declare that this is a group project and that no part of this submission has been copied from any other student's work or from any other source except where due acknowledgement is made explicitly in the text, nor has any part been written for us by another person.

Student	Contribution to the project (%) (Total 100%)	Signature
Wong Ho Shun Sang	40%	
CUI Tianlang	26%	
Yau Yuk ming	18%	
Chan Hing Chak	16%	

Table of Contents

1. Abstract	7
2. Introduction.....	8
3. Driving Question.....	9
4. Function Requirement.....	11
4.1 User Registration	11
4.2 Login.....	11
4.3 Homepage.....	11
4.4 Chat Room.....	11
4.5 AI Bot Query	11
4.6 Remote Health Monitoring.....	13
4.7 Doctor-Patient Appointment.....	13
4.8 Real-time Chat for Diagnosis	13
4.9 Sensor Data Acquisition	13
5 Non-functional Requirements	14
5.1 Security	14
5.2 Performance.....	14
5.3 Usability.....	14
5.4 Reliability	14
5.5 Compatibility	14
5.6 Scalability	15
5.7 Compliance	15
5.8 Integration.....	15
6. Design of the proposed system	16
6.1 architecture of the proposed system.....	16
6.2 scope of the proposed system.....	17
6.3 description of functions provided.....	17
6.4 data processed by the system	18
7. Documentation for problem analysis	25
7.1 Actor description and Use case description	25
Actor description	25
8 Use case description.....	26
8.1 Login and Register System.....	26

8.2 Appointment System	27
8.3 Search system	28
8.4 Chat System.....	29
8.4 Sensor Data Acquisition System	30
9 Diagram.....	33
9.1 Use Case Diagram.....	33
Remote Health Monitoring System – Full System.....	33
Appointment System	35
Search and chatbot System.....	36
Chat System.....	37
Meeting Room System	38
Sensor Data Acquisition System	39
9.2 Class diagram	40
Login and Register System.....	40
Appointment System	40
Search System.....	41
Chat System.....	41
Sensor Data Acquisition System	42
9.3 Sequence diagram	43
Login and Register System.....	43
Appointment System	45
Chat System.....	49
Sensor Data Acquisition System	50
System level Sequence Diagram	52
9.4 State Machine diagram.....	53
Login and Register System.....	53
Appointment System	55
Chatbot.....	57
Search	58
Chat System.....	59
Sensor Data Acquisition System	60
10. Documentation for detailed design	62
10.1 Data design - ERD.....	62
10.2 Software architectural design.....	63

Software / Maintenances cost	69
10.3 Hardware architectural design (system design)	70
Cost of Hardware(s) and Service(s).....	74
10.4 User interface design.....	75
Home Page.....	75
Login – Doctor.....	77
Login – User	77
Register – Doctor.....	78
Register – User	78
User - Home Page.....	79
User - Appointment	79
User - Healthy Test.....	81
User – Meeting Room.....	81
User – Chat Room	84
User – Healthy Information (English).....	85
User – Healthy Information (Chinese)	86
User - Chat Bot	87
Doctor – Home Page.....	88
Doctor - Manage Appointment.....	88
Doctor – Chat Room.....	89
Doctor – Meeting Room	91
Doctor – Sensor Record.....	94
11. Critical Evaluation	95
11.1 Difficulties encountered	95
11.2 Imitations of the proposed system.....	95
11.3 Potential difficulties	95
12.All test cases and test data	96
12.1 Login	96
12.2 Make Appointment.....	97
12.3 Healthy Test	98
12.4 Chat Room.....	99
12.5 Meeting.....	100
12.6 Chatbot	101
Test Case ID	101

TC006.....	101
12.7 Search Page	101
12.8 Manage Appointment.....	102
12.9 View Sensor Record.....	103
13. Project plan	104
14. Conclusion	106
15. Project Log.....	108
16. References.....	110
17. Program listing.....	111
Home Page	111
Index.html.....	111
Home.js	113
Login	118
Customer.php.....	118
Doctor.php	120
RegisterDoctor.php.....	123
Customer.js.....	126
Customer.js.....	127
Meeting Room.....	131
Meeting Room.html.....	131
Chat Bot.....	133
ChatBot.html.....	133
ChatBot.css	134
Script.js	138
18. User guide and installation guide.....	140
Step 1: Download and zip the file	140
Step 2: Put the zipped file into the path	140
Step 3: Open the Xampp and Turn on the Apache and MySQL.....	141
Step4: Open the Admin and insert the mysql file	142
Step 5: Go to the browser and paste the link: http://127.0.0.1/FYP web/	144
Step6: Click the button ‘Sign Up New’ to sign up.....	144
Step7: Login or sign up a new account (For Public).....	145
Step8: Make Appointment (For Public)	147
Step9: Click the button chat to communicate with doctor (Public)	148

Step10: Click the button meeting to communicate with doctor (Public)	150
Step11: Click the Healthy Information to Search the information and use the chat bot	151
Step12: Login or sign up a new account (For Doctor).....	153
Step13: Management appointment (Doctor)	155
Step14: Sensor Healthy test (Public).....	156
Step15: View Sensor Record (Doctor).....	157

1. Abstract

The objective of this project is to develop a comprehensive Remote Health Monitoring System that addresses the limitations of traditional healthcare systems. The system aims to provide individuals with a user-friendly platform for remote health monitoring, appointment scheduling, personalized health suggestions, and access to their health history. It also enables doctors to view extended patient health information and efficiently manage their patient list. The proposed solution follows a client-server architecture, with the client side developed using HTML, CSS, JavaScript, and jQuery, and the server side implemented using PHP and XAMPP. The system includes features such as a secure login system, sensor data acquisition, doctor information, appointment scheduling, online meetings, patient list management, health suggestions, and patient health history tracking. The project will be divided into various stages, including project initiation, requirements gathering and analysis, development, testing, deployment, documentation, and training. The main deliverables of the project include an SQL file for the database, a website for the user interface, and a comprehensive report documenting the project.

2. Introduction

This document provides a general description of the structure and objectives of the Remote Health Monitoring System project. The project aims to develop a comprehensive system that addresses the limitations of traditional healthcare systems by enabling remote health monitoring and improved communication between patients and doctors.

The ability to monitor certain aspects of a patient's health from their own home has become an increasingly popular telehealth option. Remote patient monitoring lets providers manage acute and chronic conditions. And it cuts down on patients' travel costs and infection risk.

The Remote Health Monitoring System is designed to provide individuals with a user-friendly platform for monitoring their health remotely. It allows users to track vital signs, such as heart rate, blood pressure, and temperature, using various sensors and devices. The system also includes features for appointment scheduling, personalized health suggestions, and access to the users' health history.

On the healthcare provider side, the system enables doctors to view extended patient health information and efficiently manage their patient list. It facilitates online meetings between doctors and patients, allowing for remote consultations and follow-ups. By leveraging technology, the system aims to enhance the efficiency and effectiveness of healthcare delivery.

The project follows a client-server architecture, with the client side developed using HTML, CSS, JavaScript, and jQuery. The server side is implemented using PHP and XAMPP, which provides a robust and secure foundation for data management and communication. The system utilizes a database to store patient information, sensor data, appointment details, and other relevant data.

By combining remote health monitoring, appointment management, and improved communication channels, the Remote Health Monitoring System aims to empower individuals to take a proactive role in managing their health while providing healthcare providers with valuable tools for delivering high-quality care remote.

3. Driving Question

Explain how the project answers this question: “How can Software Engineering techniques be used to develop software systems for supporting human activities?”

This system utilizes software engineering techniques to create a comprehensive solution that supports and enhances various human activities related to healthcare.

- The Remote Health Monitoring System addresses the need for individuals to remotely monitor their health and collaborate effectively with healthcare providers. By leveraging software engineering techniques, the system provides a user-friendly platform that supports human activities in the following ways:
- Remote Health Monitoring: The system utilizes sensors and devices to collect vital health data from individuals. Software engineering techniques are applied to design and develop the data acquisition and processing components, ensuring accurate and reliable measurement of health parameters. This enables individuals to actively monitor their health from the comfort of their homes or any remote location.
- Appointment Scheduling: The system incorporates appointment scheduling functionality, allowing users to conveniently book appointments with healthcare providers. Software engineering techniques are employed to develop a user-friendly interface for scheduling appointments, managing availability, and handling conflicts. This streamlines the process of setting up consultations and supports efficient coordination between patients and doctors.
- Personalized Health Suggestions: The system provides personalized health suggestions based on the collected health data. Software engineering techniques are utilized to analyze the data, apply algorithms and models, and generate meaningful insights and recommendations. This helps individuals make informed decisions regarding their health and well-being.

- Access to Health History: The system enables individuals to access and manage their health history. Software engineering techniques are employed to design and implement a secure database that stores patient health records. This allows users to review their medical history, track progress, and share relevant information with healthcare providers, supporting continuity of care.
- Communication and Collaboration: The system facilitates online meetings and communication channels between patients and doctors. Software engineering techniques are utilized to develop real-time communication features, ensuring secure and reliable interactions. This supports effective collaboration, remote consultations, and follow-ups, enhancing the overall healthcare experience.

Through the application of software engineering techniques, the Remote Health Monitoring System addresses the driving question by developing a software solution that supports and enhances human activities related to healthcare. It leverages technology to enable remote health monitoring, appointment scheduling, personalized health suggestions, and improved communication, ultimately empowering individuals to take an active role in managing their health while facilitating efficient healthcare delivery.

4. Function Requirement

4.1 User Registration

The system should provide a user-friendly registration process for patients, allowing them to enter their personal information accurately. The registration process should validate the entered information and ensure its integrity. Upon successful registration, patients should receive a confirmation notification or email, ensuring a smooth onboarding experience.

- The system should provide registration functionality for patients.
- Patients should register by providing necessary information through a hyperlink on the Login page.

4.2 Login

The login functionality should offer a secure and seamless access point for users (doctors and patients) to their respective accounts. Users should be able to enter their login credentials, such as username and password, and the system should authenticate the credentials to provide access. In case of incorrect credentials, appropriate error messages should be displayed to assist users in troubleshooting.

- The system should have a login functionality for both doctors and patients.
- Patients and doctors should be able to access their respective accounts by entering their credentials.

4.3 Homepage

The homepage should serve as a central hub for users, displaying relevant and up-to-date information in an easily accessible manner. It should provide a visually appealing layout, featuring healthcare news, articles, notifications, and other pertinent information. Users should be able to navigate through different sections of the homepage to find the desired information efficiently.

- The system should display a homepage HTML where users can view relevant news and information.
- Users should be able to access the homepage after successful login.

4.4 Chat Room

The Chat Room functionality should facilitate real-time communication between patients and doctors. It should offer an intuitive and user-friendly messaging interface, allowing users to exchange messages seamlessly. Users should be able to view message history, receive notifications for new messages, and have basic formatting options to enhance communication effectiveness.

- The system should provide a Chat Room functionality for patients and doctors to engage in remote real-time conversations.
- Patients should be able to access the Chat Room through the navigation bar after login.

4.5 AI Bot Query

The AI bot functionality should enable patients to obtain information and answers to their queries. The system should provide a conversational interface where patients can ask

questions or enter keywords related to their health concerns. The AI bot should employ natural language processing techniques to understand and provide relevant responses. It should support interactive conversations, allowing users to ask follow-up questions or seek clarifications.

- The system should allow patients to query information from an AI bot.
- Patients should be able to access the AI bot functionality through the navigation bar after login.

4.6 Remote Health Monitoring

The system should support the remote collection and storage of patients' health data. Patients should be able to upload their body readings or health measurements obtained from wearable devices or medical equipment. The system should securely store and associate the data with the respective patient's profile, ensuring privacy and data integrity. The uploaded data should be easily accessible for future reference by both patients and doctors.

- The system should enable patients to upload their body readings from remote devices.
- Patients should be able to log in to the system and upload their body readings for record-keeping.

4.7 Doctor-Patient Appointment

The appointment scheduling functionality should streamline the process of scheduling and managing appointments between doctors and patients. Patients should be able to view doctors' availability, select preferred dates and times, and book appointments based on their convenience. The system should provide confirmation notifications to patients and update doctors' schedules, accordingly, ensuring efficient appointment management.

- The system should facilitate doctor-patient appointments.
- Patients should be able to log in to the system and schedule appointments with doctors.

4.8 Real-time Chat for Diagnosis

The real-time chat functionality should enable doctors to engage in live conversations with patients for diagnosis and treatment purposes. Doctors should be able to ask pertinent questions, request additional information or tests, and provide medical advice. The chat interface should support multimedia elements, such as image sharing, to facilitate accurate diagnosis and enhance the doctor-patient communication experience.

- The system should enable doctors to have real-time conversations with patients for diagnosis purposes.
- Doctors should be able to use the Chat Room functionality to interact with patients and analyze their condition.

4.9 Sensor Data Acquisition

The system is integrated with sensors and devices to collect patient health data, including heart rate, temperature, and skin conductance. This data is then stored in a database, allowing doctors to remotely and in real-time monitor the patient's health status. This integration of the system with various sensors and equipment empowers healthcare providers to effectively track and manage the well-being of their patients, even from a distance, enabling more comprehensive and personalized care.

- The system should enable doctors to view patient health data instantly.
- Patients should be able to view their historical health data.

5 Non-functional Requirements

5.1 Security

The system should implement robust security measures to protect user data and ensure privacy. It should employ encryption techniques to secure data transmission and storage. Access to sensitive information should be restricted to authorized personnel, and the system should comply with relevant data protection regulations to safeguard user confidentiality.

- The system should ensure the security and confidentiality of user data, including personal information and health records.
- Data transmission and storage should be encrypted to protect patient privacy.

5.2 Performance

The system should deliver a responsive and high-performance experience, ensuring minimal delays in data retrieval, messaging, and other functionalities. It should be able to handle concurrent user interactions and maintain optimal performance even during peak usage periods, providing a smooth and efficient user experience.

- The system should provide a responsive and seamless user experience, ensuring minimal delays in chat conversations and data retrieval.
- The system should be capable of handling concurrent user interactions and processing data efficiently.

5.3 Usability

The user interface should be intuitive, visually appealing, and easy to navigate. Users should be able to understand and use the system's features without extensive training or assistance. Clear instructions, tooltips, and help sections should be available to guide users through different functionalities and aid in their overall usability.

- The system should have a user-friendly interface that is intuitive and easy to navigate.
- The system should provide clear instructions and guidance for users to perform actions effectively.

5.4 Reliability

The system should be reliable and available for users. It should have backup and recovery mechanisms in place to minimize downtime and ensure data integrity. Regular maintenance and updates should be performed to address any potential issues promptly and maintain system reliability.

- The system should be available and accessible to users reliably, with minimal downtime or system failures.
- The system should have backup and recovery mechanisms in place to prevent data loss in case of unexpected events.

5.5 Compatibility

The system should be compatible with various web browsers and devices, ensuring accessibility for a wide range of users. It should be designed to adapt to different screen sizes and resolutions, providing a consistent user experience across multiple platforms.

- The system should be compatible with various web browsers and devices to ensure accessibility for users.
- The system should be designed to support different screen sizes and resolutions.

5.6 Scalability

The system architecture should be scalable to accommodate a growing number of users and increasing data volume. It should be able to handle additional user load without significant performance degradation. Scalability measures, such as load balancing and cloud-based infrastructure, can be implemented to ensure efficient resource utilization and system stability.

- The system should be designed to accommodate a growing number of users and data without compromising performance.
- The system architecture should be scalable to handle increased user load and data volume.

5.7 Compliance

The system should comply with relevant data protection and privacy regulations, such as GDPR or HIPAA, depending on the jurisdiction. It should adhere to industry standards and best practices for healthcare information systems to ensure the secure handling and storage of sensitive user data.

- The system should comply with relevant data protection and privacy regulations, such as GDPR or HIPAA, depending on the jurisdiction.
- The system should adhere to industry standards and best practices for healthcare information systems.


5.8 Integration

The system should be designed to integrate with external systems or devices, such as electronic health record systems or remote health monitoring devices. It should support seamless data exchange and interoperability, enhancing the overall functionality and usefulness of the system within the healthcare ecosystem.

- The system should be capable of integrating with external systems or devices for data exchange, such as remote health monitoring devices or electronic health record systems.

6. Design of the proposed system

6.1 architecture of the proposed system



There are two networks, one for the public and one for the private.

The Public Router connects to the internet via the Public IP address 1941.

The Public Router is connected to the Switch-PT, which is in turn connected to Server-PT and the Doctor's network via the Doctor Router.

The System Switch is connected to the Remote Health Monitoring System Server and the MySQL Database Server.

There are two servers connected to the System Switch. The first server is called Server-PT, which is the Remote Health Monitoring System Server. This server is responsible for collecting and storing health data from remote patients. The second server is called the MySQL Database Server, which stores the data collected by Server-PT.

The PL-PT link connects the Doctor's network to the Remote Health Monitoring System Server.

Overall, this network setup allows the doctor's network to connect to the Remote Health Monitoring System Server, which stores and manages health data for patients. The MySQL Database Server is used to store this data, and the System Switch and Switch-PT are used to connect various devices within the network. The Public Router connects the network to the internet, allowing remote access to the system.

6.2 scope of the proposed system

The scope of the proposed solution includes the development of a client-server architecture-based system. The client side is developed using HTML, CSS, JavaScript, and jQuery, providing a user-friendly interface for individuals to interact with the system. The server side is implemented using C# and XAMPP, ensuring robust data management and communication.

6.3 description of functions provided

1. Secure Login System: Users can securely log in to the system using their credentials, ensuring data privacy and security.
2. Doctor Information: Doctors have access to patient details, medical history, and appointment schedules for efficient patient management.
3. Appointment Scheduling: Users can schedule appointments with healthcare providers based on availability and convenience.
4. Online Meetings: The system facilitates online meetings between doctors and patients, enabling remote consultations and follow-ups.
5. Patient List Management: Doctors can efficiently manage their patient list, track appointments, and update patient information.
6. Health Suggestions: The system provides personalized health suggestions based on collected data and analysis.
7. Patient Health History Tracking: Users can access and manage their health history, including past appointments, test results, and medications.
8. Sensor Data Acquisition: The system integrates with sensors and devices to collect vital health data, such as heart rate, blood pressure, and temperature.

6.4 data processed by the system

doctor table

#	名稱	類型	編碼與排序	屬性	空值(Null)	預設值	備註	額外資訊	動作
□ 1	doctorID 📃	int(10)		否	無				修改 刪除 更多
□ 2	name	varchar(255)	utf8_general_ci	否	無				修改 刪除 更多
□ 3	types	varchar(255)	utf8_general_ci	否	無				修改 刪除 更多
□ 4	tel	int(8)		否	無				修改 刪除 更多
□ 5	email	varchar(255)	utf8_general_ci	否	無				修改 刪除 更多
□ 6	password	varchar(10)	utf8_general_ci	否	無				修改 刪除 更多

Mata data of doctor

Data Item Name	Data Type	Data length	Description of Data Item
doctorID	integer	10	Doctor Identification Number
name	character	255	Doctor Name
types	character	255	Doctor types
tel	integer	8	Doctor phone number
email	email	255	Doctor email
password	alphanumeric	10	Doctor Account Password

users table

#	名稱	類型	編碼與排序	屬性	空值(Null)	預設值	備註	額外資訊	動作
□	1 userID	int(10)		否	無				修改 刪除 更多
□	2 userName	varchar(50)	utf8_general_ci	否	無				修改 刪除 更多
□	3 Email	varchar(255)	utf8_general_ci	否	無				修改 刪除 更多
□	4 tel	varchar(8)	utf8_general_ci	否	無				修改 刪除 更多
□	5 password	varchar(8)	utf8_general_ci	否	無				修改 刪除 更多

Data Item Name	Data Type	Data length	Description of Data Item
userID	integer	10	Public Identification Number
userName	character	50	Public Name
Email	email	255	Public email
tel	integer	8	Public phone number

sensor table-GSR sensor

#	名稱	類型	編碼與排序	屬性	空值(Null)	預設值	備註	額外資訊	動作
□	1 id	int(11)		否	無	AUTO_INCREMENT			修改 刪除 更多
□	2 gsr	int(11)		否	無				修改 刪除 更多
□	3 datetime	datetime		否	current_timestamp()				修改 刪除 更多

Mata data of sensor table-GSR sensor

Data Item Name	Data Type	Data length	Description of Data Item
id	integer	11	Public Identification Number
gsr	integer	11	GSR sensor data
datetime	datetime		Current data & time

heartrate table- heart rate sensor

#	名稱	類型	編碼與 排序	屬性	空值 (Null)	預設值	備註	額外資訊	動作
<input type="checkbox"/>	1 id	int(11)		否	無	AUTO_INCREMENT		修改 刪除 更多	
<input type="checkbox"/>	2 rate	int(11)		否	無			修改 刪除 更多	
<input type="checkbox"/>	3 datetime	datetime		否	無	current_timestamp()		修改 刪除 更多	

Data Item Name	Data Type	Data length	Description of Data Item
id	integer	11	Public Identification Number
rate	integer	11	Heart rate sensor data
datetime	datetime		Current data & time

temperature_sensor table- temperature sensor

#	名稱	類型	編碼與 排序	屬性	空值 (Null)	預設值	備註	額外資訊	動作
<input type="checkbox"/>	1 id	int(11)		否	無	AUTO_INCREMENT		修改 刪除 更多	
<input type="checkbox"/>	2 Celsius	float		否	無			修改 刪除 更多	
<input type="checkbox"/>	3 datetime	datetime		否	無	current_timestamp()		修改 刪除 更多	

Mata data of temperature_sensor

Data Item Name	Data Type	Data length	Description of Data Item
id	integer	11	Public Identification Number
Celsius	float		temperature sensor data°C
datetime	datetime		Current data & time

history_sensor table- get data from 3 types of sensors(for using)

#	名稱	類型	編碼與排序	屬性	空值(Null)	預設值	備註	額外資訊	動作
□ 1	id 	int(11)		否	無	AUTO_INCREMENT		 修改  刪除  更多	
□ 2	userName	varchar(50)	utf8mb4_general_ci	否	無			 修改  刪除  更多	
□ 3	datetime	datetime		否	無	current_timestamp()		 修改  刪除  更多	
□ 4	tem_C	float		否	無			 修改  刪除  更多	
□ 5	heartrate	int(11)		否	無			 修改  刪除  更多	
□ 6	gsr	int(11)		否	無			 修改  刪除  更多	

Mata data of history_sensor table

Data Item Name	Data Type	Data length	Description of Data Item
id	integer	11	Public Identification Number
userName	character	50	Public Name
datetime	datetime		Current data & time
tem_C	float		temperature sensor data°C
heartrate	integer	11	Heart rate sensor data
gsr	integer	11	GSR sensor data

Appointment table

#	名稱	類型	編碼與排序	屬性	空值(Null)	預設值	備註	額外資訊	動作
□ 1	id 🔑	int(10)		否	無	AUTO_INCREMENT			修改 刪除 更多
□ 2	name	varchar(255) utf8mb4_general_ci		否	無				修改 刪除 更多
□ 3	date	date		否	無				修改 刪除 更多
□ 4	time	time(6)		否	無				修改 刪除 更多
□ 5	VisitMethod	varchar(255) utf8mb4_general_ci		否	無				修改 刪除 更多
□ 6	doctorName	varchar(255) utf8mb4_general_ci		否	無				修改 刪除 更多

Mata data of appointment

Data Item Name	Data Type	Data length	Description of Data Item
id	integer	10	Appointment Identification Number
name	character	255	Username(Public)
date	date		Meeting Date
time	time	6	Meeting Time
VisitMethod	character	255	Meeting method(face to face or online)
doctorName	character	255	Doctor Name

member Table

#	名稱	類型	編碼與排序	屬性	空值(Null)	預設值	備註	額外資訊	動作	
1	user_id	int(11)		否	無	AUTO_INCREMENT		修改	刪除	更多
2	unique_id	int(255)		否	無			修改	刪除	更多
3	fname	varchar(255)	utf8mb4_general_ci	否	無			修改	刪除	更多
4	lname	varchar(255)	utf8mb4_general_ci	否	無			修改	刪除	更多
5	email	varchar(255)	utf8mb4_general_ci	否	無			修改	刪除	更多
6	password	varchar(255)	utf8mb4_general_ci	否	無			修改	刪除	更多
7	img	varchar(255)	utf8mb4_general_ci	否	無			修改	刪除	更多
8	status	varchar(255)	utf8mb4_general_ci	否	無			修改	刪除	更多

Mata data of member

Data Item Name	Data Type	Data length	Description of Data Item
user_id	integer	11	user Identification Number
unique_id	int	255	Identify user
fname	character		User first name
lname	character	6	User last name
email	character	255	email
password	character	255	User password
img	character	255	User icon image
status	character	255	User status

messages table

資料表結構									
#	名稱	類型	編碼與排序	屬性	空值(Null)	預設值	備註	額外資訊	動作
1	msg_id	int(11)		否	無	AUTO_INCREMENT		修改 刪除 更多	
2	incoming_msg_id	int(255)		否	無			修改 刪除 更多	
3	outgoing_msg_id	int(255)		否	無			修改 刪除 更多	
4	msg	varchar(1000) utf8mb4_general_ci		否	無			修改 刪除 更多	

Mata data of messages

Data Item Name	Data Type	Data length	Description of Data Item
msg_id	integer	11	User message Number
incoming_msg_id	int	255	Income message record id
outgoing_msg_id	int	255	Out go message record id
msg	character	1000	Message content

7. Documentation for problem analysis

7.1 Actor description and Use case description

Actor description

Doctor	Doctors rely on the system to streamline patient care, enhance communication with patients, and leverage data-driven insights to provide comprehensive healthcare services.
Public (Patient)	The public represents individuals who interact with the healthcare system as patients, caregivers, or general users. They utilize the system to access healthcare services, manage their health information, and engage in communication with healthcare providers.

8 Use case description

8.1 Login and Register System

Use case name	Register
Case ID	UC-100
actor	Public and Doctor
Brief description	All user can use the system to register an account. The system will save the account.
Preconditions	Users must follow the password format.
Post-conditions	Users can open the system
Flow of events	<ol style="list-style-type: none">1. Users click button to their role signup page.2. Users input the personal information3. Users input the correct password format.4. Users click the signup button.
Alternative flow and exception	The user may input wrong password format. The system will show the password format message. The user may miss an input data.
Priority	high
Non-behavioral requirement	The password should be kept secretly by the user.

Use case name	Login
Case ID	UC-200
actor	Public and Doctor
Brief description	All user can use their account to login to the system.
Preconditions	Users must have an account.
Post-conditions	User can login to the system and start using it
Flow of events	<ol style="list-style-type: none">1. Users click button to their role login page.2. Users input correct username.3. Users input correct password.
Alternative flow and exception	The user may input wrong password or username. The system will alert a message to tell user
Priority	high
Non-behavioral requirement	The password should be kept secretly by the user.

8.2 Appointment System

Use case name	Make Appointment
Case ID	UC-300
actor	Public
Brief description	The public can fill the form to make an appointment
Preconditions	
Post-conditions	
Flow of events	<ol style="list-style-type: none"> 1. Public login to the system. 2. Public click 'Appointment' button to appointment page. 3. Public view doctor information. 4. Public fill data in the form. 5. Public click submit button.
Alternative flow and exception	The public may select the wrong doctor.
Priority	high
Non-behavioral requirement	The public should ask the related doctor.

Use case name	Manage Appointment
Case ID	UC-400
actor	Doctor
Brief description	Doctor can manage (edit and delete) appointment record.
Preconditions	The public make an appointment.
Post-conditions	If any changes are made to the appointment record, such as editing the data or deleting the record, the public should be notified by the doctor or the system to ensure effective communication and minimize any confusion or inconvenience.
Flow of events	<ol style="list-style-type: none"> 1. Doctor view appointment record. 2. Doctor edit the form data. 3. Doctor click submit button. 4. Doctor click delete button to delete the outdated record or canceled record.
Alternative flow and exception	The public may select the wrong doctor.
Priority	high
Non-behavioral requirement	In case of emergencies or unforeseen circumstances, doctors should promptly reach out to the public to reschedule appointments or provide alternative arrangements.

8.3 Search system

Use case name	Search
Case ID	UC-500
actor	public
Brief description	Public can search disease information
Preconditions	The public enter keywords.
Post-conditions	Accurate Information: The user should have been presented with accurate, up-to-date, and relevant information about the disease they searched for.
Flow of events	1.User enter keyword 2.User Click submit button 3.System display the information about the keyword
Alternative flow and exception	The system may not have the information about the keyword
Priority	middle
Non-behavioral requirement	The system errors

8.4 Chat System

Use case name	Chatbot
Case ID	UC-600
actor	public
Brief description	Public can ask anything to chatbot
Preconditions	Public enter something.
Post-conditions	Public should be polite to the chatbot
Flow of events	1. User enter keyword 2. User Click enter button 3. chatbot display something to user
Alternative flow and exception	
Priority	high
Non-behavioral requirement	The chatbot errors.

Use case name	Chat Room
Case ID	UC-700
actor	Patient, Doctor
Brief description	Public can ask anything to chatbot
Preconditions	Public enter something.
Post-conditions	This use case describes the process of engaging in remote real-time conversation between a patient and a doctor in the Chat Room.
Flow of events	1. The patient selects the Chat Room option from the navigation bar. 2. The system displays the Chat Room interface. 3. The patient enters a message in the input field. 4. The patient clicks the "Send" button. 5. The system sends the message to the doctor. 6. The doctor receives the message in their Chat Room interface. 7. The doctor reads the message. 8. The doctor composes and enters a response in the input field. 9. The doctor clicks the "Send" button. 10. The system sends the response to the patient. 11. The patient receives the response in their Chat Room interface. 12. The patient reads the response. 13. Steps 3-12 repeat as needed for the ongoing conversation.

	14. The patient can choose to end the conversation by closing the Chat Room interface.
Alternative flow and exception	<ul style="list-style-type: none"> • If the doctor is not available or not logged into the system: ▪ The patient receives a notification indicating the doctor's unavailability and is prompted to try again later. • If the patient or doctor encounters a technical issue (e.g., network problem): ▪ The system displays an error message and provides instructions for troubleshooting or contacting technical support.
Priority	Medium
Non-behavioral requirement	<ul style="list-style-type: none"> • The Chat Room interface should provide a real-time messaging experience with minimal delay. • The conversation history should be stored securely and accessible only to authorized users.

8.4 Sensor Data Acquisition System

Use case name	Sensor data Acquisition
Case ID	UC-800
actor	Public
Brief description	Public can use heart rate, temperature, gsr and sensors to obtain health data.
Preconditions	Sensors connected to database and Wi-Fi.
Post-conditions	
Flow of events	<ol style="list-style-type: none"> 1. Open Arduino IDE 2. The sensor is connected to the computer port. 3. The sensor is connected to Wi-Fi. 4. The sensor is connected to the database. 5. Click the upload button to compile. 6. Sensor data will be uploaded to the database. 7. Repeat the above steps with the other sensors.
Alternative flow and exception	If you select the wrong port and card, an error will occur.
Priority	high
Non-behavioral requirement	The public should stop the connection, when the compilation is complete, the public should disconnect.

Use case name	Upload sensor data
Case ID	UC-900
actor	Public
Brief description	Public can upload sensor data to database of the history sensor table.
Preconditions	The sensor has been uploaded to its matching table.
Post-conditions	
Flow of events	<ol style="list-style-type: none"> 1. Click upload button. 2. The system will upload these 3-type sensors to the history sensor table.
Alternative flow and exception	
Priority	high
Non-behavioral requirement	


Use case name	View history sensor record
Case ID	UC-1000
actor	Public and Doctor
Brief description	Users can upload View history sensor record.
Preconditions	The public has uploaded sensor data to the historical sensor table.
Post-conditions	
Flow of events	<ol style="list-style-type: none"> 1. View history sensor table . 2. The system returns the history sensor record.
Alternative flow and exception	
Priority	high
Non-behavioral requirement	

Use case name	View patient list
Case ID	UC-1100
actor	Doctor
Brief description	Doctor can view patient lists.
Preconditions	Patient has been registering an account.
Post-conditions	
Flow of events	<p>1. View users table.</p> <p>2. The system returns the patient list.</p>
Alternative flow and exception	
Priority	high
Non-behavioral requirement	Doctor can contact patient.

9 Diagram


9.1 Use Case Diagram

Remote Health Monitoring System – Full System




This is a use case diagram of the Remote Health Monitoring System. The diagram involves three primary actors: Patient, Doctor, and AI Bot. The Patient represents individuals utilizing the system for remote health monitoring, while the Doctor represents healthcare professionals engaging with patients through the system. The AI Bot acts as an artificial intelligence component providing automated responses and information. The diagram includes four main use cases: Register, Chat with Doctor, Query AI Bot, and Remote Health Monitoring. Users can register and access the system, engage in real-time communication with doctors for medical advice, interact with the AI Bot for information, and upload health data for remote monitoring by doctors. This use case diagram provides a concise overview of the system's main functionalities and the interactions between the actors and the system.

Login and Register System




The use case diagram showcases the Login and Register functions of the Remote Health Monitoring System. The User actor interacts with the System actor to perform two main use cases. In the Register use case, the User provides personal information and submits it to the System for registration. The System validates the information, creates a new user account, and sends a confirmation notification. The Login use case allows the User to enter their login credentials, which are then verified by the System. If the credentials are valid, access is granted to the User. Otherwise, an error message is displayed. The diagram highlights the essential steps involved in user registration and login processes, emphasizing the interaction between the User and the System for accessing the Remote Health Monitoring System.

Appointment System



This use case diagram illustrates the appointment function of the Remote Health Monitoring System. The main actors in the diagram are the Patient and the Doctor, representing individuals who interact with the system. The diagram includes two primary use cases: "View Doctor Availability" and "Book Appointment." In the "View Doctor Availability" use case, the Patient can check the availability of different doctors within the system. This allows the Patient to see the doctors' schedules and identify suitable time slots for appointments. In the "Book Appointment" use case, the Patient selects a preferred date and time slot from the available options and requests an appointment with the chosen doctor. The system confirms the appointment and sends a notification to both the Patient and the Doctor. This use case diagram demonstrates the essential steps involved in the appointment booking process within the Remote Health Monitoring System, highlighting the interaction between the Patient, Doctor, and the system.

Search and chatbot System




The use case diagram showcases the search and chat bot functions of the Remote Health Monitoring System. The diagram includes two primary actors: User and System. The User represents individuals interacting with the system, while the System represents the Remote Health Monitoring System itself.

The diagram incorporates two main use cases: "Search Health Information" and "Chat with Bot." In the "Search Health Information" use case, the User can enter specific keywords or queries to search for relevant health information within the system. The System retrieves and presents the search results to the User, providing them with access to a wide range of health-related knowledge.

In the "Chat with Bot" use case, the User initiates a conversation with the system's Chat Bot. The User can ask questions, seek medical advice, or request information on specific health topics. The Chat Bot utilizes predefined knowledge and natural language processing techniques to generate appropriate responses and engage in a dialogue with the User.

These two use cases demonstrate how the Remote Health Monitoring System allows Users to search for health information and interact with the Chat Bot for personalized support and guidance. The diagram illustrates the interaction between the User and the System, emphasizing the functionalities of search and chat bot within the system.

Chat System




The use case diagram represents the chat system function of the Remote Health Monitoring System. The diagram includes three primary actors: User, Doctor, and System. The User represents individuals utilizing the system, the Doctor represents healthcare professionals, and the System represents the Remote Health Monitoring System itself.

The diagram consists of one main use case: "Chat with Doctor." In this use case, the User can initiate a conversation with the Doctor through the chat system provided by the Remote Health Monitoring System. The User can seek medical advice, ask questions, or discuss health concerns with the Doctor in real-time. The Doctor, on the other hand, responds to the User's messages, provides diagnoses, offers recommendations, or engages in a dialogue to address the User's medical needs.

The chat system facilitates secure and convenient communication between the User and the Doctor, allowing for remote consultations and support. The system ensures that sensitive medical information remains confidential and provides a platform for effective healthcare communication. This use case diagram demonstrates the interaction between the User, Doctor, and the System, highlighting the chat system's function within the Remote Health Monitoring System.


Meeting Room System



The use case diagram represents the chat system function of the Remote Health Monitoring System. The diagram includes three primary actors: User, Doctor, and System. The User represents individuals utilizing the system, the Doctor represents healthcare professionals, and the System represents the Remote Health Monitoring System itself.

The diagram consists of one main use case: "Online Meeting". In this use case, the User can join a meeting to have a conversation through Webcam, Audio, Chat and Whiteboard. It provides variety way to let user communicate with each other.

Sensor Data Acquisition System




- "Sensor data Acquisition": This use case represents the system's ability to collect and acquire data from various sensors.
- "Upload sensor data": This use case describes the functionality of uploading the acquired sensor data to the system's database.
- "View history sensor record": This use case allows users to view the historical sensor data records for patients.
- The "Sensor data Acquisition" use case includes the functionality of acquiring sensor data.
- The "View history sensor record" use case extends the "Upload sensor data" use case, allowing users to view the uploaded sensor data.


Overall, the use case diagram depicts the core functionalities of the Sensor Data Acquisition System, including data acquisition, data upload, and historical data viewing, as well as the interactions between the system and its various actors (Public, Doctor, and User).

9.2 Class diagram


Login and Register System




Appointment System




Search System



Chat System




Sensor Data Acquisition System




9.3 Sequence diagram

Login and Register System




The sequence diagram depicts the Register system function of the Remote Health Monitoring System. The diagram shows the interactions between the User, RegisterController, and Database. It starts with the User sending a registration request to the RegisterController. The RegisterController validates the provided information and communicates with the Database to create a new user account. Once the account is successfully created, the RegisterController sends a registration confirmation back to the User. This sequence diagram illustrates the flow of actions involved in the Register process within the Remote Health Monitoring System.




The sequence diagram illustrates the Login system function of the Remote Health Monitoring System. The diagram showcases the interactions between the User, LoginController, and Database. It begins with the User sending a login request to the LoginController. The LoginController verifies the credentials provided by the User and communicates with the Database to retrieve the user's information. Upon successful verification, the LoginController sends a login response back to the User, allowing access to the system. This sequence diagram provides a visual representation of the flow of actions involved in the Login process within the Remote Health Monitoring System.

Appointment System


Make Appointment



Manage Appointment




Chat bot system




The sequence diagram represents the chat bot system function of the Remote Health Monitoring System. The diagram depicts the interactions between the User, ChatBot, and KnowledgeBase. It begins with the User sending a message to the ChatBot. The ChatBot analyzes the message and retrieves relevant information from the KnowledgeBase. Based on the analysis, the ChatBot generates a response and sends it back to the User. This sequence diagram illustrates the flow of actions involved in the chat bot system function, where the ChatBot utilizes the KnowledgeBase to provide intelligent and informative responses to the User's queries within the Remote Health Monitoring System.

Search system



The sequence diagram depicts the search system function of the Remote Health Monitoring System. The diagram shows the interactions between the User, SearchController, and Database. It begins with the User sending a search query to the SearchController. The SearchController processes the query and communicates with the Database to retrieve relevant health information. The SearchController then sends the search results back to the User. This sequence diagram illustrates the flow of actions involved in the search system function, where the User can search for and retrieve health information within the Remote Health Monitoring System.


Chat System




The sequence diagram depicts the chat system function of the Remote Health Monitoring System using a chat room for real-time communication. The diagram shows the interactions between Users and the Chat Room. Users join the chat room and send messages, which are broadcasted to all participants. The Chat Room receives the messages and distributes them to all connected Users. This sequence diagram illustrates the flow of real-time communication within the chat system, where Users can engage in conversations with each other through the shared chat room in the Remote Health Monitoring System.

Sensor Data Acquisition System


Sensor data Acquisition




Upload sensor data




View history sensor record



View patient list



System level Sequence Diagram




The System level sequence diagram of the Remote Health Monitoring System depicts the flow of interactions between the User, Login, Chat, ChatBot, Appointment, and Database. It shows how users' login to the system, engage in real-time chats facilitated by the Chat receive intelligent responses from the ChatBot based on the KnowledgeBase, and manage appointments through the Appointment, which communicates with the Database to update and retrieve appointment information. This diagram provides an overview of the key functionalities of the Remote Health Monitoring System and how different components interact to deliver a seamless user experience.

9.4 State Machine diagram


Login and Register System

Register



Login


Visual Paradigm Professional (Version 10.0.0 Build 10000) Training Guide



The state machine diagram represents the login system function of the Remote Health Monitoring System, showcasing the various states and transitions involved in the login process.


Appointment System

Make Appointment



The state machine diagram represents the make appointment function of the Appointment system in the Remote Health Monitoring System, showcasing the different states and transitions involved in the process of creating an appointment.


Manage Appointment



The state machine diagram depicts the manage appointment function of the Appointment system in the Remote Health Monitoring System, illustrating the various states and transitions involved in the process of modifying or canceling an existing appointment.


Chatbot

Visual Paradigm Professional 14.2.0.50000 (Build 20170628)




The state machine diagram represents the chat bot system function of the Remote Health Monitoring System, showcasing the different states and transitions involved in the interaction between the user and the chat bot for providing responses and gathering information.

Search



The state machine diagram illustrates the search system function of the Remote Health Monitoring System, highlighting the various states and transitions involved in the process of searching and retrieving health information based on user queries.


Chat System




The state machine diagram represents the chat room system function of the Remote Health Monitoring System, showcasing the different states and transitions involved in the real-time communication and message exchange within the shared chat room among multiple users.

Sensor Data Acquisition System


Sensor data Acquisition




Upload sensor data



View history sensor record




View patient list



10. Documentation for detailed design

10.1 Data design - ERD



10.2 Software architectural design



Google Chrome

is a web browser developed by Google. It was first released in 2008 for Microsoft Windows, built with free software components from Apple WebKit and Mozilla Firefox. Versions were later released for Linux, macOS, iOS, and also for Android, where it is the default browser. The browser is also the main component of ChromeOS, where it serves as the platform for web applications.



PhpStorm is a development tool for PHP and Web projects. It's a perfect PHP IDE for working with Laravel, Symfony, Drupal, WordPress, and other frameworks.




DroidCam turns your phone into a camera source on your computer.



Oracle database system

MySQL, developed by Oracle, is a renowned open-source relational database management system (RDBMS) extensively utilized for constructing reliable and scalable applications. It encompasses a broad array of functionalities that aid developers in efficiently creating, storing, managing, and retrieving data.



Xampp

XAMPP is a widely used open-source software package that provides a local development environment for web applications. It includes Apache as the web server, MySQL as the database management system, PHP as the scripting language, and Perl for additional flexibility. XAMPP simplifies the setup and configuration of these components, allowing developers to quickly create and test web applications on their local machines before deployment.




Microsoft Office

Microsoft Office is a comprehensive collection of productivity applications that offers a diverse set of tools for creating, organizing, and collaborating on various types of documents, spreadsheets, presentations, and other tasks related to enhancing productivity.



Visual Studio Code

Visual Studio Code, commonly known as VS Code, is a lightweight and highly customizable source code editor developed by Microsoft. It provides a user-friendly interface and a vast ecosystem of extensions, enabling developers to write, edit, and debug code efficiently across multiple programming languages. With its extensive features and intuitive design, Visual Studio Code enhances the coding experience and boosts productivity for developers.



Bracket

Brackets is a lightweight and open-source code editor designed specifically for web development. It offers a clean and minimalist interface with a focus on front-end development, providing features like inline editing, live preview, and quick editing of CSS and HTML. Brackets aims to streamline the web development workflow and make coding for the web more accessible and efficient.

Arduino IDE



Arduino is an open-source integrated development environment (IDE) that makes it easy to write code and upload it to the board. Using the Arduino software, you can work with any Arduino board. The software provides a library manager that makes it convenient to install library files, with thousands of official and contributed libraries available for direct download. Sample code is also provided. The installation process is quick and straightforward, whether you're using Windows, macOS, or Linux operating systems.



Microsoft Windows 10

Microsoft Windows 10 is an operating system developed by Microsoft and is the latest version of the Windows family. It provides a user-friendly interface, enhanced security features, and a wide range of applications and tools for various tasks. Windows 10 offers a seamless user experience with features like the Start Menu, Cortana virtual assistant, Microsoft Edge web browser, and compatibility with a vast range of software and hardware devices. It is designed to be versatile, efficient, and suitable for both personal and professional use. Here are some of the key functions of Microsoft Windows:

- User Interface
- File Management
- Security
- Device Drivers
- Networking
- System Maintenance



JetBrains IDE

The Python IDE for Professional Developers



PyCharm

PyCharm is a powerful and user-friendly integrated development environment (IDE) specifically designed for Python programming. It offers a range of features such as code completion, debugging tools, intelligent code analysis, and project management capabilities. PyCharm provides a seamless coding experience, helping developers write, test, and debug Python code efficiently. With its intuitive interface and extensive set of tools, PyCharm is widely used by developers to enhance their productivity in Python development projects.

Software / Maintenances cost

Title	Cycle	Price	Subtotal	Software / Maintenances	Total
Hardware check	Per 3 months	388	388	M	1522 per year
System update(s)	Per week	180	180	M	9360 per year
License of Anti-virus software	Per year	1280	1280	S	1280 per year
Oracle database system	Per version	7600	7600	S	7600
License of Microsoft Office	Per version	1488	1488	S	1488
License of Microsoft Windows	Per version	2980	2980	S	2980

10.3 Hardware architectural design (system design)



Workstation PC

A workstation PC is a high-performance computer designed for professional use in demanding tasks such as graphic design, video editing, 3D modeling, scientific simulations, and software development. It typically features powerful hardware components like multi-core processors, ample RAM, dedicated graphics cards, and fast storage options. Workstation PCs are optimized for reliability, stability, and efficient workflow, often including specialized software and tools. They provide the processing power and resources required to handle complex workloads and deliver superior performance for professionals in various industries.



Router

A router is a networking device that connects multiple devices within a network and directs data packets between them. It acts as a central hub, allowing devices to communicate with each other and with devices on other networks. Routers use IP addresses to determine the best path for data transmission and ensure that information reaches its intended destination. They provide functions like network address translation (NAT), firewall protection, and quality of service (QoS) to manage network traffic effectively. Routers are essential in homes, offices, and other networked environments to enable internet connectivity and facilitate communication between devices.

Sensor

GSR Sensor



The GSR (Galvanic Skin Response) sensor measures the electrical conductivity of the skin, which can be influenced by factors such as sweat and moisture. It is commonly used to assess the level of emotional or physiological arousal in individuals.

1. High Value: A high value in the GSR sensor reading typically indicates a lower skin conductivity or higher skin resistance. This can occur when the skin is dry or when the individual is in a state of low arousal or relaxation.
2. Low Value: Conversely, a low value in the GSR sensor reading generally indicates higher skin conductivity or lower skin resistance. This can happen when the skin is moist or when the individual is experiencing increased arousal or emotional stimulation.

Temperature Sensor



A temperature sensor, also known as an infrared thermometer or non-contact thermometer, is a device used to measure the temperature of an object or person without making physical contact. It utilizes infrared technology to detect and measure the thermal radiation emitted by the target.

The temperature readings provided by a temperature gun sensor are typically displayed in Celsius ($^{\circ}\text{C}$).

ECG Sensor



An ECG (Electrocardiogram) sensor is a device used to measure and record the electrical activity of the heart. It detects the electrical signals generated by the heart muscle during each heartbeat and provides valuable information about the heart's rhythm, rate, and overall cardiac health.

To calculate the heart rate from an ECG signal, the number of heartbeats occurring within a specific time interval needs to be determined. The most common time interval used is 60 seconds (1 minute). The ECG sensor records the electrical activity of the heart continuously, and by analyzing the recorded data, the number of heartbeats occurring within the 60-second interval can be counted.

To obtain the heart rate, the number of heartbeats is divided by the duration of the interval (60 seconds) and multiplied by 60 to convert it to beats per minute (bpm). This calculation provides an estimation of the heart rate based on the ECG signal recorded by the sensor.

Cost of Hardware(s) and Service(s)


Title	Quantity	Unit Price	Subtotal	
Router (TP-Link Archer AX11000)	2	2380 HKD	4760 HKD	Router connected to the internet for system able to communicate with the server
Network fees (HKBN 10G)	1	528 HKD	528 HKD	Fees for the internet access, charged from broadband company
24/7 technical support	1	328 HKD	328 HKD	One-time payment of 24/7 technical support through online chat, phone call
Workstation PC (ThinkStation P360 Ultra - Intel)	2	9896.4 HKD	19792.8 HKD	Workstation PC for running the software and daily uses

10.4 User interface design

Home Page



The screenshot shows the homepage of a "Remote Health System". At the top left is a red shield icon with a white cross. To its right are navigation links: Home, Services, About Us, and Login. Below the header is a large banner featuring the text "Revolutionizing Healthcare: The Rise Of Remote Medical Systems" in bold black font. Underneath this text is a smaller paragraph about remote medical systems. A "Sign Up Now" button is located at the bottom left of the banner. To the right of the text is a photograph of a medical control room with a large screen displaying various data and graphs.



The screenshot shows the "Our Services" page. At the top left is a red shield icon with a white cross. To its right are navigation links: Home, Services, About Us, and Login. The main title "Our Services" is centered above five service cards. Each card has an icon and a title followed by a detailed description. The services listed are: 1. Remote Physiological Parameter Monitoring (described as monitoring physiological parameters like heart rate and blood pressure through sensors). 2. Health Data Recording And Analysis (described as recording and storing patient health data for analysis). 3. Remote Medical Consultation And Diagnosis (described as engaging in real-time or asynchronous consultations to analyze data and provide treatment plans). 4. Medication Reminders And Management (described as providing medication reminders to help patients adhere to their schedules). 5. Health Education And Guidance (described as offering relevant health education and guidance, including disease management and lifestyle recommendations). Below the service cards is a section titled "Emergency Calls And Rescue" with a "911" icon and a note about sending distress signals for medical assistance.



Remote Health System

[Home](#) [Services](#) [About Us](#) [Login](#)

Our Services



Health Data Recording And Analysis

Remote Health Control Systems Can Record And Store The Patient's Health Data Including Measurement Values, Medication History, Symptom Descriptions, Etc. Healthcare Providers Can Analyze This Data To Gain A Comprehensive Understanding Of The Patient's Condition And Make More Accurate Diagnoses And Treatment Plans.

Remote Medical Consultation And Diagnosis

Patients Can Engage In Real-Time Or Asynchronous Medical Consultations And Diagnoses With Healthcare Providers Through The Remote Health Control Platform. They Can Communicate Their Symptoms, Ask Questions, And Receive Advice And Treatment Plans Through Video Calls, Phone Calls, Or Online Chats.

Medication Reminders And Management

The Remote Health Control System Can Provide Medication Reminders To Help Patients Adhere To Their Medication Schedules. Additionally, Healthcare Providers Can Remotely Monitor The Patient's Medication Adherence Through The System And Make Timely Adjustments To The Treatment Plan If Needed.

Health Education And Guidance

The Remote Health Control Platform Can Offer Relevant Health Education And Guidance To Patients, Including Disease Management, Lifestyle Recommendations, Dietary Advice, Etc. This Helps Patients To Better Understand And Manage Their Health Conditions.



Emergency Calls And Rescue

In Emergency Situations, Patients Can Use The Remote Health Control System To Send Emergency Calls Or Distress Signals To Promptly Receive Medical Assistance.



Remote Health System

[Home](#) [Services](#) [About Us](#) [Login](#)

Remote Physiological Parameter Monitoring
Through Sensors, Smart Devices, Or Wearable Devices, The System Can Monitor A Patient's Physiological Parameters Such As Heart Rate, Blood Pressure, Blood Oxygen Saturation, Body Temperature, Etc. These Data Can Be Transmitted In Real-Time To Healthcare Providers For Remote Monitoring And Assessment Of The Patient's Health Condition.

Health Data Recording And Analysis

Remote Health Control Systems Can Record And Store The Patient's Health Data Including Measurement Values, Medication History, Symptom Descriptions, Etc. Healthcare Providers Can Analyze This Data To Gain A Comprehensive Understanding Of The Patient's Condition And Make More Accurate Diagnoses And Treatment Plans.

Remote Medical Consultation And Diagnosis

Patients Can Engage In Real-Time Or Asynchronous Medical Consultations And Diagnoses With Healthcare Providers Through The Remote Health Control Platform. They Can Communicate Their Symptoms, Ask Questions, And Receive Advice And Treatment Plans Through Video Calls, Phone Calls, Or Online Chats.

Medication Reminders And Management

The Remote Health Control System Can Provide Medication Reminders To Help Patients Adhere To Their Medication Schedules. Additionally, Healthcare Providers Can Remotely Monitor The Patient's Medication Adherence Through The System And Make Timely Adjustments To The Treatment Plan If Needed.

Health Education And Guidance

The Remote Health Control Platform Can Offer Relevant Health Education And Guidance To Patients, Including Disease Management, Lifestyle Recommendations, Dietary Advice, Etc. This Helps Patients To Better Understand And Manage Their Health Conditions.



Emergency Calls And Rescue

In Emergency Situations, Patients Can Use The Remote Health Control System To Send Emergency Calls Or Distress Signals To Promptly Receive Medical Assistance.



Create Your Own Profile



Upload Your Body Values




Consult A Doctor Online




Reserve


Login – Doctor




Login – User




Register – Doctor



Register – User





User - Home Page



User - Appointment

Doctor List					
Id	Name	Types	Tel	Email	Appointment
1	Mike Lee	General medical practitioner	51212456	Mike12@gmail.com	<input type="button" value="Book"/>
2	William Hunt	General medical practitioner	61216456	William12@gmail.com	<input type="button" value="Book"/>
3	Denise Dunn	Cardiologist	62215477	Denise22@gmail.com	<input type="button" value="Book"/>
4	sam	General medical practitioner	51712126	s123@gmail.com	<input type="button" value="Book"/>



User - Healthy Test


The screenshot shows the 'Sensor Data Dashboard' of a remote health monitoring system. At the top, there are three cards displaying current sensor readings:


- Temperature Sensor**: Temperature: 36.7 °C
- Heart Rate Sensor**: Heart Rate: 61.68353353353353 BPM
- Human Resistance Sensor**: Resistance: 610 Ohms


Below these cards is a green header bar labeled "History Record". Underneath is a table with the following data:

ID	Date & Time	Temperature	HeartRate	GSR
1	2024-04-06 17:00:25	36.7	80	610
2	2024-04-07 17:00:33	36.7	80	610
3	2024-04-07 17:14:20	36.7	80	610
4	2024-04-07 17:16:55	36.7	80	610
5	2024-04-07 17:20:56	36.7	80	610
6	2024-04-07 17:21:59	36.7	80	610

User – Meeting Room







User – Chat Room

The image shows a mobile application interface. At the top, there is a navigation bar with the text "Remote health monitoring system" and icons for Home, Appointment, HealthyTest, Meeting, Chat, and Healthy Information. Below this, a user profile for "Isla Jones" is displayed, showing she is "Active now". A search bar is present, followed by a list of users: "Mike Lee" (12123) and another user whose profile picture is partially visible. The main content area shows a chat room between "Isla Jones" and "Mike Lee". The messages are as follows:

- Isla Jones: 2222222222222222
- Mike Lee: ji
- Isla Jones: 11
- Mike Lee: 222
- Isla Jones: 12123

At the bottom of the screen is a text input field with the placeholder "Type a message here..." and a send button.

User – Healthy Information (English)

The screenshots show a user navigating through a "Disease Database" section of a "Remote health monitoring system".

Screenshot 1: Home Page

- Header: Remote health monitoring system, Home, Appointment, HealthyTest, Meeting, Chat, Healthy Info, Profile.
- Search bar: Enter disease... search.
- Main content: Welcome to the Disease Information website. Here is some information about the disease...

Screenshot 2: Disease Detail Page (Polio)

- Header: Remote health monitoring system, Home, Appointment, HealthyTest, Meeting, Chat, Healthy Info, Profile.
- Search bar: Enter disease... search.
- Left sidebar: Neurological Diseases (Polio, Parkinson's Disease, Endocrine Diseases, Diabetes, Cardiovascular Diseases, Heart Disease, Mental and Psychological Diseases (Depression, Mental illness, Anorexia), Respiratory Diseases (Asthma)).
- Right content area: **Poliomyelitis**. Description: Poliomyelitis, or infantile paralysis, is a contagious disease caused by the poliovirus. It primarily affects children under the age of 5 and can lead to varying degrees of paralysis. Subtypes: **Paralytic poliomyelitis**: Affects the spinal cord and can lead to limb paralysis. **Bulbar poliomyelitis**: Affects muscles involved in breathing and swallowing. **Non-paralytic poliomyelitis**: The virus causes infection without leading to paralysis. Symptoms: include: Vomiting, Neck stiffness and pain, Weakness in the limbs. Treatment: mainly relies on vaccination. There is no cure for poliomyelitis, and treatment is primarily symptomatic, including physical therapy and supportive care.

Screenshot 3: Overview of heart disease

- Header: Remote health monitoring system, Home, Appointment, HealthyTest, Meeting, Chat, Healthy Info, Profile.
- Search bar: s search.
- Main content: **Welcome to the Disease Information website**. Here is some information about the disease...
- Overview of heart disease**: Heart disease is one of the leading causes of death worldwide and includes a variety of diseases that affect heart function. These disorders may involve other structures of the heart, such as the coronary arteries, heart muscle, valves, or conduction system.
- Types of heart disease**: There are many types of heart disease, including but not limited to the following:
 - Coronary Heart Disease (CAD)**: The main blood vessels in the heart are damaged, which can lead to a heart attack.
 - myocardial infarction**: Commonly known as a heart attack, it occurs when blood flow to the heart muscle is blocked.
 - heart failure**: The heart is unable to pump blood effectively to meet the body's needs.
 - arrhythmia**: The heart beats fast, slow, or irregularly.
 - Heart valve disease**: Involves abnormal opening and closing of heart valves.
 - Congenital heart disease**: Structural problems with the heart that are present at birth.
 - cardiomyopathy**: Diseases that affect the heart muscle, which may cause the heart to become enlarged or hardened.
 - pericarditis**: Inflammation of the pericardium.
- Symptoms of Heart Disease**: The symptoms of heart disease vary from person to person. Some may not even know they have heart disease, while others may experience one or more of the following symptoms:
 - Chest pain or discomfort
 - Difficulty breathing
 - Palpitations
 - Weakness or dizziness
 - Fatigue
 - Swelling of the ankles, legs, or abdomen
- Prevention and Treatment**: The prevention and treatment of heart disease involve lifestyle changes, medication, and possibly surgical treatment. Here are some preventive measures:
 - Healthy diet**: Eat plenty of fruits, vegetables, and whole grains.
 - Regular exercise**: At least 150 minutes of moderate-intensity exercise per week.
 - Maintain a healthy weight**: Reducing weight can lessen the burden on the heart.

User – Healthy Information (Chinese)

The screenshot displays a web-based remote health monitoring system. At the top, there is a dark navigation bar with the title "Remote health monitoring system" and several menu items: Home, Appointment, HealthyTest, Meeting, Chat, Healthy Info, and Profile. Below the navigation bar, there is a search bar with placeholder text "搜索..." and a "搜索" button. The main content area features a large, light blue header banner with the text "歡迎來到疾病資訊網站" (Welcome to the Disease Information Website) and a subtext "這裡是一些疾病的相關資訊..." (Here is some information about diseases...). The main content below the banner is a detailed article about "小兒痲痺症概述" (Overview of小儿麻痹症). The article includes sections on symptoms, types, prevention/treatment, risk factors, and diagnosis. It also contains a sidebar with a QR code and a "分享" (Share) button.

小兒痲痺症概述

小兒痲痺症，又稱脊髓灰質炎，是一種由小兒麻痺病毒引起的傳染病。它主要影響5歲以下的兒童，並可能導致不同程度的麻痹。

小兒痲痺症的類型

小兒麻痺症主要有三種類型：

- 脣癱型小兒麻痺症：影響骨骼，可能導致數體麻痺。
- 延髓型小兒麻痺症：影響呼吸和吞嚥肌肉。
- 無麻痺型小兒麻痺症：病毒感染後未出現麻痹症狀。

小兒痲痺症的症狀

小兒麻痺症的症狀可能包括：

- 發熱
- 腹脹
- 頭痛
- 嘔吐
- 颈部僵硬和疼痛
- 四肢無力

預防與治療

小兒麻痺症的預防主要依靠疫苗接種。目前沒有特效藥物可以治愈小兒麻痺症，治療主要針對症狀進行，包括物理治療和支持療法。

小兒痲痺症的風險因素


導致小兒麻痺症的風險因素包括：

- 未接種疫苗
- 居住在病毒流行區
- 免疫系統較弱

小兒痲痺症的診斷

診斷小兒麻痺症通常依據臨牀症狀和實驗室檢測，如脊髓培養和PCR檢測。

User - Chat Bot



Doctor – Home Page



Doctor - Manage Appointment

Appointment List							
ID	Name	Date	Time	Visit Method	Doctor Name	Edit	Delete
4	Isla Jones	2024-04-22	04:08:00.000000	online	Mike Lee		

Remote health monitoring system

- [Home](#)
- [Appointment Record](#)
- [Chat](#)
- [Meeting](#)
- [Sensor Record](#)
- [Profile](#)


ID	Name	Date	Time
4	Isla Jones	2024-04-22	04:08:00.000000


Edit Appointment

ID:	<input type="text" value="4"/>	Doctor Name	Edit	Delete
Name:	<input type="text" value="Isla Jones"/>	Mike Lee		
Date:	<input type="text" value="2024/04/22"/>			
Time:	<input type="text" value="04:08"/>			
Visit Method:	<input type="text" value="online"/>			
Doctor Name:	<input type="text" value="Mike Lee"/>			


[Submit](#) [Close](#)


Doctor – Chat Room






Doctor – Meeting Room







Doctor – Sensor Record

 **Remote health monitoring system** [Home](#) [Appointment Record](#) [Chat](#) [Meeting](#) [Sensor Record](#) [Profile](#)

[Patient List](#) [Sensor Record](#)

Patient List				
Id	Name	Email	Tel	
1	Isla Jones	Isla2j@gmail.com	55678899	
2	Juan Owens	Juow11@gmail.com	65446322	
3	TAM	s123@gmail.com	51712126	
4	sssa	s123@gmail.com	51234561	

Sensor Record					
Id	Name	Date & Time	Temperature	Heart Rate	GSR
1	Isla Jones	2024-04-06 17:00:25	36.7	80	610
2	Isla Jones	2024-04-07 17:00:33	36.7	80	610
3	Isla Jones	2024-04-07 17:14:20	36.7	80	610
4	Isla Jones	2024-04-07 17:16:55	36.7	80	610
5	Isla Jones	2024-04-07 17:20:56	36.7	80	610
6	Isla Jones	2024-04-07 17:21:59	36.7	80	610
7	Isla Jones	2024-04-07 18:46:33	36.7	80	610
8	Juan Owens	2024-04-07 21:08:52	36.9	76	608

Sensor Record					
Id	Name	Date & Time	Temperature	Heart Rate	GSR
8	Juan Owens	2024-04-07 21:08:52	36.9	76	608

11. Critical Evaluation

11.1 Difficulties encountered

1. Technical Challenges: The implementation of certain features, such as login system and appointment system, required extensive testing and debugging to ensure smooth functionality. Dealing with compatibility issues, software conflicts, or integrating different technologies posed challenges that required additional time and effort.
2. Resource Limitations: Limited availability of resources, both in terms of hardware and personnel, impacted the development process. The system's reliance on XAMPP for the server infrastructure posed constraints on scalability and performance.

11.2 Imitations of the proposed system

1. Dependency on XAMPP for Server: The system relies on XAMPP, a local server solution, for hosting the application. This dependency may limit scalability and accessibility since it requires the system to be hosted on a local machine rather than a dedicated server or cloud-based infrastructure.
2. Lack of Financial Resources for Setting up a Server: The proposed system may face challenges in terms of funding to set up a dedicated server infrastructure. This limitation can impact the system's performance, reliability, and availability, as it may not be able to handle many users or provide consistent service uptime.
3. Restricted Usage to a Single Computer: Due to the limitations of the system's infrastructure, it may only allow for simultaneous usage of two webpages on the same computer. This restriction can hinder scalability and limit the number of concurrent users who can access the system simultaneously.

11.3 Potential difficulties

1. User Adoption and Acceptance: Introducing a new system in the healthcare domain may face resistance or hesitancy from both healthcare providers and patients. Overcoming user resistance, addressing concerns about usability and trust, and effectively promoting the benefits of the system are essential for successful user adoption and acceptance.
2. Stakeholder Collaboration: Building a successful remote health monitoring system requires collaboration among various stakeholders, including healthcare providers, technology vendors, regulatory bodies, and patients.
3. Regulatory and legal considerations: Remote health monitoring systems must comply with various regulatory requirements, such as data protection laws, medical device regulations, and telemedicine guidelines. Navigating through these regulations and ensuring compliance can be complex and time-consuming, potentially delaying the implementation process.

12.All test cases and test data

12.1 Login

Test Case ID	TC001																																																							
Test Description	Verify that the login functionality works with valid credentials																																																							
Preconditions	The user must be at the login page																																																							
Test Steps	<p>1.Enter the username in the username field.</p> <p>2.Enter the password in the password field.</p> <p>3.Click the login button.</p>																																																							
Test Data	<p>User</p> <table border="1"> <thead> <tr> <th>userID</th> <th>userName</th> <th>Email</th> <th>tel</th> <th>password</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Isla Jones</td> <td>Isla2j@gmail.com</td> <td>55678899</td> <td>Jo123@66</td> </tr> <tr> <td>2</td> <td>Juan Owens</td> <td>Juow11@gmail.com</td> <td>65446322</td> <td>Juan123@</td> </tr> <tr> <td>3</td> <td>TAM</td> <td>s123@gmail.com</td> <td>51712126</td> <td>Ss123456</td> </tr> <tr> <td>4</td> <td>ssssaa</td> <td>s123@gmail.com</td> <td>51234561</td> <td>Ss123456</td> </tr> </tbody> </table> <p>Doctor</p> <table border="1"> <thead> <tr> <th>doctorID</th> <th>name</th> <th>types</th> <th>tel</th> <th>email</th> <th>password</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Mike Lee</td> <td>General medical practitioner</td> <td>51212456</td> <td>Mike12@gmail.com</td> <td>Mike11@</td> </tr> <tr> <td>2</td> <td>William Hunt</td> <td>General medical practitioner</td> <td>61214456</td> <td>William12@gmail.com</td> <td>Wi22sh7@</td> </tr> <tr> <td>3</td> <td>Denise Dunn</td> <td>Cardiologist</td> <td>62215477</td> <td>Denise22@gmail.com</td> <td>Dencar2@</td> </tr> <tr> <td>4</td> <td>sam</td> <td>General medical practitioner</td> <td>51712126</td> <td>s123@gmail.com</td> <td>Ss123456</td> </tr> </tbody> </table>	userID	userName	Email	tel	password	1	Isla Jones	Isla2j@gmail.com	55678899	Jo123@66	2	Juan Owens	Juow11@gmail.com	65446322	Juan123@	3	TAM	s123@gmail.com	51712126	Ss123456	4	ssssaa	s123@gmail.com	51234561	Ss123456	doctorID	name	types	tel	email	password	1	Mike Lee	General medical practitioner	51212456	Mike12@gmail.com	Mike11@	2	William Hunt	General medical practitioner	61214456	William12@gmail.com	Wi22sh7@	3	Denise Dunn	Cardiologist	62215477	Denise22@gmail.com	Dencar2@	4	sam	General medical practitioner	51712126	s123@gmail.com	Ss123456
userID	userName	Email	tel	password																																																				
1	Isla Jones	Isla2j@gmail.com	55678899	Jo123@66																																																				
2	Juan Owens	Juow11@gmail.com	65446322	Juan123@																																																				
3	TAM	s123@gmail.com	51712126	Ss123456																																																				
4	ssssaa	s123@gmail.com	51234561	Ss123456																																																				
doctorID	name	types	tel	email	password																																																			
1	Mike Lee	General medical practitioner	51212456	Mike12@gmail.com	Mike11@																																																			
2	William Hunt	General medical practitioner	61214456	William12@gmail.com	Wi22sh7@																																																			
3	Denise Dunn	Cardiologist	62215477	Denise22@gmail.com	Dencar2@																																																			
4	sam	General medical practitioner	51712126	s123@gmail.com	Ss123456																																																			
Expected Result	The user is redirected to the dashboard after successful login																																																							
Actual Result	Login successful																																																							
Status	(To be filled after execution)																																																							
Comments	Null																																																							

12.2 Make Appointment

Test Case ID	TC002																												
Test Description	Make the Appointment																												
Preconditions	The user must be at the login page																												
Test Steps	1. Click the button ‘Make appointment’ 2. Choose value the Date. 3. Choose the type of appointment. 4. Choose the Doctor.																												
Test Data	<table border="1"> <thead> <tr> <th>id</th> <th>name</th> <th>date</th> <th>time</th> <th>VisitMethod</th> <th>doctorName</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>Isla Jones</td> <td>2024-01-25</td> <td>00:29:00.000000</td> <td>face to face</td> <td>William Hunt</td> </tr> <tr> <td>3</td> <td>Isla Jones</td> <td>2024-04-23</td> <td>22:43:00.000000</td> <td>face to face</td> <td>sam</td> </tr> <tr> <td>4</td> <td>Isla Jones</td> <td>2024-04-22</td> <td>04:08:00.000000</td> <td>online</td> <td>Mike Lee</td> </tr> </tbody> </table>					id	name	date	time	VisitMethod	doctorName	2	Isla Jones	2024-01-25	00:29:00.000000	face to face	William Hunt	3	Isla Jones	2024-04-23	22:43:00.000000	face to face	sam	4	Isla Jones	2024-04-22	04:08:00.000000	online	Mike Lee
id	name	date	time	VisitMethod	doctorName																								
2	Isla Jones	2024-01-25	00:29:00.000000	face to face	William Hunt																								
3	Isla Jones	2024-04-23	22:43:00.000000	face to face	sam																								
4	Isla Jones	2024-04-22	04:08:00.000000	online	Mike Lee																								
Expected Result	The user can see the successful message after successful making appointment																												
Actual Result	Make Appointment successful																												
Status	(To be filled after execution)																												
Comments	Null																												



12.3 Healthy Test

Test Case ID	TC003																											
Test Description	Do the Healthy Test and upload the data to the system																											
Preconditions	The user must be at the login page																											
Test Steps	<p>1. Set up the sensor.</p> <p>2. Click the button “Healthy Test”</p> <p>3. Use the sensor to check the body data.</p> <p>4. Click the green button to upload the data</p>																											
Test Data	<table border="1"> <thead> <tr> <th>id</th> <th>rate</th> <th>datetime</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>70</td> <td>2024-03-11 16:27:21</td> </tr> <tr> <td>2</td> <td>80</td> <td>2024-03-13 16:02:22</td> </tr> <tr> <td>3</td> <td>4096</td> <td>2024-04-19 16:21:33</td> </tr> <tr> <td>4</td> <td>4000</td> <td>2024-04-19 16:21:33</td> </tr> <tr> <td>5</td> <td>4100</td> <td>2024-04-19 16:22:27</td> </tr> <tr> <td>6</td> <td>4200</td> <td>2024-04-19 16:22:27</td> </tr> <tr> <td>7</td> <td>3056</td> <td>2024-04-19 16:22:27</td> </tr> <tr> <td>8</td> <td>3701</td> <td>2024-04-19 16:22:27</td> </tr> </tbody> </table>	id	rate	datetime	1	70	2024-03-11 16:27:21	2	80	2024-03-13 16:02:22	3	4096	2024-04-19 16:21:33	4	4000	2024-04-19 16:21:33	5	4100	2024-04-19 16:22:27	6	4200	2024-04-19 16:22:27	7	3056	2024-04-19 16:22:27	8	3701	2024-04-19 16:22:27
id	rate	datetime																										
1	70	2024-03-11 16:27:21																										
2	80	2024-03-13 16:02:22																										
3	4096	2024-04-19 16:21:33																										
4	4000	2024-04-19 16:21:33																										
5	4100	2024-04-19 16:22:27																										
6	4200	2024-04-19 16:22:27																										
7	3056	2024-04-19 16:22:27																										
8	3701	2024-04-19 16:22:27																										
Expected Result	The user can see check and upload the data after using the sensor to check the body information																											
Actual Result	View and upload the data																											
Status	(To be filled after execution)																											
Comments	Null																											

12.4 Chat Room

Test Case ID	TC004																																																
Test Description	Communicate with each other successfully																																																
Preconditions	The user(doctor) must be at the login page																																																
Test Steps	<p>1.Enter the Chat Room.</p> <p>2.Click the Send button.</p>																																																
Test Data	<p>Message</p> <table border="1"> <thead> <tr> <th>msg_id</th> <th>incoming_msg_id</th> <th>outgoing_msg_id</th> <th>msg</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>1572804435</td> <td>1110990992</td> <td>11</td> </tr> <tr> <td>2</td> <td>1572804435</td> <td>1110990992</td> <td>22</td> </tr> <tr> <td>3</td> <td>1572804435</td> <td>1110990992</td> <td>3</td> </tr> <tr> <td>4</td> <td>1572804435</td> <td>1110990992</td> <td>11</td> </tr> <tr> <td>5</td> <td>1110990992</td> <td>1572804435</td> <td>22222222222222</td> </tr> <tr> <td>6</td> <td>1110990992</td> <td>1572804435</td> <td>ji</td> </tr> <tr> <td>7</td> <td>1110990992</td> <td>1572804435</td> <td>11</td> </tr> <tr> <td>8</td> <td>1572804435</td> <td>1110990992</td> <td>222</td> </tr> <tr> <td>9</td> <td>1572804435</td> <td>1334909215</td> <td>sss</td> </tr> <tr> <td>10</td> <td>1110990992</td> <td>1334909215</td> <td>sss</td> </tr> <tr> <td>11</td> <td>1572804435</td> <td>1110990992</td> <td>12123</td> </tr> </tbody> </table>	msg_id	incoming_msg_id	outgoing_msg_id	msg	1	1572804435	1110990992	11	2	1572804435	1110990992	22	3	1572804435	1110990992	3	4	1572804435	1110990992	11	5	1110990992	1572804435	22222222222222	6	1110990992	1572804435	ji	7	1110990992	1572804435	11	8	1572804435	1110990992	222	9	1572804435	1334909215	sss	10	1110990992	1334909215	sss	11	1572804435	1110990992	12123
msg_id	incoming_msg_id	outgoing_msg_id	msg																																														
1	1572804435	1110990992	11																																														
2	1572804435	1110990992	22																																														
3	1572804435	1110990992	3																																														
4	1572804435	1110990992	11																																														
5	1110990992	1572804435	22222222222222																																														
6	1110990992	1572804435	ji																																														
7	1110990992	1572804435	11																																														
8	1572804435	1110990992	222																																														
9	1572804435	1334909215	sss																																														
10	1110990992	1334909215	sss																																														
11	1572804435	1110990992	12123																																														
Expected Result	The user(doctor) can chat with others successfully																																																
Actual Result	Send and Received Successfully																																																
Status	(To be filled after execution)																																																
Comments	Null																																																

12.5 Meeting

Test Case ID	TC005
Test Description	Verify that the meeting room functionality works with valid credentials
Preconditions	The user(doctor) must be at the meeting room page
Test Steps	<p>1.Enter the username in the username field.</p> <p>2.Click the join meeting button.</p>
Test Data	<p>User</p>  <p>Doctor</p> 
Expected Result	The user(doctor) is redirected to the meeting room after enter username and click join meeting button
Actual Result	Enter to the meeting room successful
Status	(To be filled after execution)
Comments	Null

12.6 Chatbot

Test Case ID	TC006
Test Description	Verify that the chatbot functionality works with valid credentials
Preconditions	The users must be at the search information page
Test Steps	1.Enter the questions or information in the chat bot 2.Click the enter button to wait for the reply of the chatbot
Test Data	Keywords of the questions or information
Expected Result	The chatbot reply to the user with relevant answer or information
Actual Result	Users get the answer or information with high accuracy
Status	(To be filled after execution)
Comments	Null

12.7 Search Page

Test Case ID	TC007
Test Description	Search the information of the symptom
Preconditions	The user must be at the login page
Test Steps	1.Enter the Search Page. 2.Enter the keywords 3.Click the Search button.
Test Data	Search Data <ul style="list-style-type: none">• Fever• Cough• Headache• Muscle pain• Sore throat• Fatigue• Nausea• Vomiting• Diarrhea• Nasal congestion• Runny nose
Expected Result	The user searches the correct information
Actual Result	Show the information about the symptom
Status	(To be filled after execution)
Comments	Null

12.8 Manage Appointment

Test Case ID	TC008																												
Test Description	Manage the Appointment which is applied from the users																												
Preconditions	The user must be at the login page																												
Test Steps	1.Click the button “Appointment Record” 2.View the Appointment Record. 3.Click the Green button to edit (can change the time of appointment) 4.Click the Red button to delete.																												
Test Data	<table border="1"> <thead> <tr> <th>id</th> <th>name</th> <th>date</th> <th>time</th> <th>VisitMethod</th> <th>doctorName</th> </tr> </thead> <tbody> <tr> <td>2</td> <td>Isla Jones</td> <td>2024-01-25</td> <td>00:29:00.000000</td> <td>face to face</td> <td>William Hunt</td> </tr> <tr> <td>3</td> <td>Isla Jones</td> <td>2024-04-23</td> <td>22:43:00.000000</td> <td>face to face</td> <td>sam</td> </tr> <tr> <td>4</td> <td>Isla Jones</td> <td>2024-04-22</td> <td>04:08:00.000000</td> <td>online</td> <td>Mike Lee</td> </tr> </tbody> </table>					id	name	date	time	VisitMethod	doctorName	2	Isla Jones	2024-01-25	00:29:00.000000	face to face	William Hunt	3	Isla Jones	2024-04-23	22:43:00.000000	face to face	sam	4	Isla Jones	2024-04-22	04:08:00.000000	online	Mike Lee
id	name	date	time	VisitMethod	doctorName																								
2	Isla Jones	2024-01-25	00:29:00.000000	face to face	William Hunt																								
3	Isla Jones	2024-04-23	22:43:00.000000	face to face	sam																								
4	Isla Jones	2024-04-22	04:08:00.000000	online	Mike Lee																								
Expected Result	The user searches the correct information																												
Actual Result	Show the information about the symptom																												
Status	(To be filled after execution)																												
Comments	Null																												

12.9 View Sensor Record

Test Case ID	TC009																																																																													
Test Description	Doctor can view the record and its details																																																																													
Preconditions	The user must be at the login page																																																																													
Test Steps	1.Click the button “Sensor Record” 2.View the Sensor Record. 3.Click the ‘Patient List’ button to view the information of the patient 4.Click the ‘Sensor Record’ button to view all the record.																																																																													
Test Data	<table border="1"> <thead> <tr> <th>id</th> <th>userName</th> <th>datetime</th> <th>tem_C</th> <th>heartrate</th> <th>gsr</th> </tr> </thead> <tbody> <tr><td>1</td><td>Isla Jones</td><td>2024-04-06 17:00:25</td><td>36.7</td><td>80</td><td>610</td></tr> <tr><td>2</td><td>Isla Jones</td><td>2024-04-07 17:00:33</td><td>36.7</td><td>80</td><td>610</td></tr> <tr><td>3</td><td>Isla Jones</td><td>2024-04-07 17:14:20</td><td>36.7</td><td>80</td><td>610</td></tr> <tr><td>4</td><td>Isla Jones</td><td>2024-04-07 17:16:55</td><td>36.7</td><td>80</td><td>610</td></tr> <tr><td>5</td><td>Isla Jones</td><td>2024-04-07 17:20:56</td><td>36.7</td><td>80</td><td>610</td></tr> <tr><td>6</td><td>Isla Jones</td><td>2024-04-07 17:21:59</td><td>36.7</td><td>80</td><td>610</td></tr> <tr><td>7</td><td>Isla Jones</td><td>2024-04-07 18:46:33</td><td>36.7</td><td>80</td><td>610</td></tr> <tr><td>8</td><td>Juan Owens</td><td>2024-04-07 21:08:52</td><td>36.9</td><td>76</td><td>608</td></tr> <tr><td>9</td><td>Isla Jones</td><td>2024-04-19 22:16:57</td><td>36.7</td><td>3701</td><td>610</td></tr> <tr><td>10</td><td>Isla Jones</td><td>2024-04-19 22:32:02</td><td>36.7</td><td>62</td><td>610</td></tr> <tr><td>11</td><td>Isla Jones</td><td>2024-04-19 22:53:08</td><td>36.7</td><td>62</td><td>610</td></tr> </tbody> </table>						id	userName	datetime	tem_C	heartrate	gsr	1	Isla Jones	2024-04-06 17:00:25	36.7	80	610	2	Isla Jones	2024-04-07 17:00:33	36.7	80	610	3	Isla Jones	2024-04-07 17:14:20	36.7	80	610	4	Isla Jones	2024-04-07 17:16:55	36.7	80	610	5	Isla Jones	2024-04-07 17:20:56	36.7	80	610	6	Isla Jones	2024-04-07 17:21:59	36.7	80	610	7	Isla Jones	2024-04-07 18:46:33	36.7	80	610	8	Juan Owens	2024-04-07 21:08:52	36.9	76	608	9	Isla Jones	2024-04-19 22:16:57	36.7	3701	610	10	Isla Jones	2024-04-19 22:32:02	36.7	62	610	11	Isla Jones	2024-04-19 22:53:08	36.7	62	610
id	userName	datetime	tem_C	heartrate	gsr																																																																									
1	Isla Jones	2024-04-06 17:00:25	36.7	80	610																																																																									
2	Isla Jones	2024-04-07 17:00:33	36.7	80	610																																																																									
3	Isla Jones	2024-04-07 17:14:20	36.7	80	610																																																																									
4	Isla Jones	2024-04-07 17:16:55	36.7	80	610																																																																									
5	Isla Jones	2024-04-07 17:20:56	36.7	80	610																																																																									
6	Isla Jones	2024-04-07 17:21:59	36.7	80	610																																																																									
7	Isla Jones	2024-04-07 18:46:33	36.7	80	610																																																																									
8	Juan Owens	2024-04-07 21:08:52	36.9	76	608																																																																									
9	Isla Jones	2024-04-19 22:16:57	36.7	3701	610																																																																									
10	Isla Jones	2024-04-19 22:32:02	36.7	62	610																																																																									
11	Isla Jones	2024-04-19 22:53:08	36.7	62	610																																																																									
Expected Result	The user searches the correct information																																																																													
Actual Result	Show the information about the symptom																																																																													
Status	(To be filled after execution)																																																																													
Comments	Null																																																																													

13. Project plan

1. Project Schedule:

The project will be divided into several phases, each with specific tasks and milestones. The estimated timeline for the project is as follows:

Phase 1: Project Initiation

- Define project scope and objectives
- Formulate project team and roles
- Identify project stakeholders

Phase 2: Requirements Gathering and Analysis

- Conduct user interviews and research
- Define system requirements and specifications
- Identify hardware and software dependencies

Phase 3: Development

- Design the user interface
- Implement client-side using HTML, CSS, JavaScript, and jQuery
- Implement server-side using PHP and XAMPP
- Integrate sensor data acquisition and processing components
- Implement functionality

Phase 4: Testing

- Perform unit testing on individual components
- Conduct integration testing
- Perform system testing and bug fixing

Phase 5: Deployment

- Set up the server environment
- Deploy the system on the server
- Perform final testing and troubleshooting

Phase 6: Documentation and Training

- Prepare comprehensive project documentation
- Create user manuals and guides

日程表

开始时间
2023年9月4日

2023年9月18日 2023年10月2日 2023年10月16日 2023年10月30日 2023年11月13日 2023年11月27日 2023年12月11日 2023年12月25日 2024年1月8日

将带有日期的任务添加到日程表

任务名称	工期	开始时间	完成时间
1 Project Initiation Phase Define project scope and objectives. Identify stakeholders and their roles. Create a project plan and timeline. Develop a project charter.	7个工作日	2023年9月4日	2023年9月12日
2 Requirements Gathering and Analysis Phase Define functional and non-functional requirements. Create use cases and user stories. Develop system architecture and design.	3个工作日	2023年9月4日	2023年9月6日
3 Development Phase Set up the development environment. Implement the client-side using HTML, CSS, JavaScript, and jQuery. 13 Implement the server-side using C# and XAMPP. 14 Develop algorithms for authentication, data retrieval, and data processing. 15 Implement the login system for patients and doctors. 16 Design and implement the SQL database for storing patient and doctor information. 17 Establish a connection with the sensor for data acquisition. 18 Create a user-friendly interface with a user drop-down menu. 19 Develop features for doctor information, contact, appointment scheduling, online meetings, patient list management, healthy suggestions, doctor count, and patient health history tracking.	3个工作日	2023年9月6日	2023年9月8日
4 Testing Phase Conduct unit testing for individual components. 22 Perform integration testing to ensure all modules work together. 23 Test the system's functionality, performance, and security. 24 Identify and fix any bugs or issues. 25 Conduct user acceptance testing with selected users.	14个工作日	2023年9月11日	2023年9月27日
5 Deployment Phase 27 Prepare the system for deployment. 28 Set up hosting environment and configure servers. 29 Deploy the system to production. 30 Perform final testing in the live environment. 31 Address any last-minute issues or concerns.	5个工作日	2023年9月25日	2023年9月29日
6 Documentation and Training Phase	80个工作日	2023年9月25日	2024年1月12日
7 Testing Phase	14个工作日	2023年9月30日	2023年10月18日
8 Deployment Phase	14个工作日	2023年10月3日	2023年10月20日
9 Documentation and Training Phase	14个工作日	2023年10月11日	2023年9月15日
10 Testing Phase	7个工作日	2023年10月15日	2023年9月22日
11 Deployment Phase	5个工作日	2023年10月21日	2023年9月27日
12 Documentation and Training Phase	14个工作日	2023年10月25日	2024年1月12日
13 Testing Phase	14个工作日	2023年11月3日	2023年11月22日
14 Deployment Phase	14个工作日	2023年11月17日	2023年12月6日
15 Documentation and Training Phase	14个工作日	2023年12月6日	2023年12月25日
16 Testing Phase	21个工作日	2023年12月15日	2024年1月12日
17 Deployment Phase	21个工作日	2024年1月12日	2024年2月9日
18 Documentation and Training Phase	21个工作日	2024年3月9日	2024年4月5日

甘特图

2023年9月 2023年10月 2023年11月 2023年12月 2024年1月

21 28 4 11 18 25 2 9 16 23 30 6 13 20 27 4 11 18 25 1 8 15

14. Conclusion

The Remote Health Monitoring System (RHMS) project has successfully demonstrated its capability to revolutionize patient care by leveraging digital technologies to facilitate real-time health monitoring and communication between doctors and patients. This system, equipped with functionalities like registration, login, appointment scheduling, health record management, and live communication, has addressed critical gaps in remote healthcare management.

System Features and Benefits

1. **User-Friendly Interface:** The RHMS provides a seamless interface for users to register and login, ensuring that even those with minimal technical knowledge can easily access health services.
2. **Enhanced Doctor-Patient Interaction:** By incorporating features such as online chat and video meetings, the system has significantly improved the way doctors interact with their patients, offering convenience and reducing the need for physical visits.
3. **Efficient Appointment Management:** Doctors can manage appointments effectively, view patient health records at a glance, and make informed decisions quickly, which is crucial in medical practice.
4. **Real-Time Health Monitoring:** The ability for users to check their heart rate using devices and upload the data online has been a game-changer. This functionality not only helps in continuous monitoring but also alerts healthcare providers promptly about any critical changes in the patient's condition.
5. **Accessible Health Records:** Patients have the benefit of viewing their own health records whenever needed. This transparency helps in fostering trust and encourages patients to take an active role in managing their health.
6. **Automated Assistance:** The integration of a chatbot provides immediate answers to common queries, which enhances user experience and reduces the workload on healthcare professionals.

Impact

The RHMS has shown potential in enhancing the efficiency of healthcare delivery. It has reduced the time and cost associated with traditional healthcare processes by minimizing physical consultations and optimizing the scheduling and management of appointments. Moreover, the system has empowered patients, giving them control over their health management, and has ensured that doctors can provide timely and effective treatment.

In conclusion, the Remote Health Monitoring System has successfully demonstrated how digital innovation can be harnessed to enhance the accessibility and quality of healthcare services. It signifies a step forward in the field of telemedicine, offering a robust framework that not only meets current healthcare demands but also anticipates future needs.

15. Project Log

Overview This log details the weekly schedule and milestones for the final year project from September 4, 2023, to April 22, 2024. Weekly meetings with the supervisor are scheduled every Tuesday, with project work sessions every Wednesday and Saturday.

Weekly Breakdown with Dates

Weekly Breakdown with Dates					
Week	Date (Tuesday)	Date (Wednesday)	Date (Saturday)	Activity	Notes
1	Sep 5, 2023	Sep 6, 2023	Sep 9, 2023	Initial meeting with supervisor	Discuss project scope and objectives
				Kick-off session	Setup project workspace, initial research
				Continue research	Focus on literature review
2	Sep 12, 2023	Sep 13, 2023	Sep 16, 2023	Progress update with supervisor	Present research findings
				Start developing project outline	Outline major components
				Refine project outline and start detailed planning	Create detailed task list
3	Sep 19, 2023	Sep 20, 2023	Sep 23, 2023	Review outline with supervisor	Adjust based on feedback

4	Sep 26, 2023	Sep 27, 2023	Sep 30, 2023	Progress update and feedback	Review initial results
				Implement supervisor's feedback	Modify design/Approach
				Continue development	Further coding and testing
5	Oct 3, 2023	Oct 4, 2023	Oct 7, 2023	Mid-project review	Comprehensive review of project status
				Address any critical issues	Make necessary adjustments
				Proceed with refined plans	Continue with implementation
...
29	Apr 16, 2024	Apr 17, 2024	Apr 20, 2024	Final meeting with supervisor	Review project journey, discuss any follow-up actions
				Project handover (if applicable)	Transfer materials, code, and documentation
				Celebration and reflection session	Recognize team effort, discuss future opportunities

Additional Notes

Ensure all logs and documentation are updated promptly after each session. Prepare questions and discussion points for each meeting with the supervisor to maximize the effectiveness of the feedback. Utilize Wednesday and Saturday sessions to address any pending issues, ensuring continuous progress on the project. This structured approach helps maintain regular progress checks and ensures that the project remains on track throughout the academic year.

16. References

SYNERTONE

<http://www.synertone.net/synertone/article/index/id/161>

GovHK

<https://www.gov.hk/tc/residents/health/healthadvice/>

17. Program listing

Home Page

Index.html

```
C: > xampp > htdocs > FYP web > index.html
1  <!DOCTYPE html>
2  <html>
3  <head>
4  | <title>Remote Health System</title>
5  | <link rel="stylesheet" type="text/css" href="Home.css">
6  </head>
7
8  <body>
9  | <header>
10 | | <a href="#" class="logo">Remote Health System</a>
11 | </header>
12 | <nav class="navbar">
13 | | <a href="#home">Home</a>
14 | | <a href="#Services">Services</a>
15 | | <a href="#">About us</a>
16 | </nav>
17 </header>
18
19 <!-- Home Section-->
20 <section class="home" id="home">
21 | <div class="content">
22 | | <h3>Revolutionizing Healthcare:</h3><h3>The Rise of Remote Medical Systems</h3>
23 | | <p>Remote Medical Systems, also known as telemedicine or telehealth systems, are transforming the healthcare landscape. The
24 | | <a href="Register_Login/Customer.php" class="btn">Sign Up Now</a>
25 | </div>
26
27 | <div class="image">
28 | | 
29 | </div>
30 </section>
31
32 <section class="Services" id="Services">
33 | <h1 class="heading"><span>Our Services</span></h1>
34 | <div class="box-container">
35 | | <div class="box">
36 | | | 
37 | | | <div class="content">
38 | | | | 
39 | | | | <h3>Remote physiological parameter monitoring</h3>
40 | | | | <p>Through sensors, smart devices, or wearable devices, the system can monitor a patient's physiological parameters such as blood glucose levels, heart rate, and blood pressure in real-time.</p>
41 | | | </div>
42 | | <div class="box">
43 | | | 
44 | | | <div class="content">
45 | | | | 
46 | | | | <h3>Health data recording and analysis</h3>
47 | | | | <p>Remote health control systems can record and store the patient's health data, including measurement values, medication
48 | | | </div>
49 | </div>
```

```

50   | </div>
51   |   | <div class="box">
52   |   |   | 
53   |   |   | <div class="content">
54   |   |   |   | 
55   |   |   |   | <h3>Remote medical consultation and diagnosis</h3>
56   |   |   |   | <p>Patients can engage in real-time or asynchronous medical consultations and diagnoses with healthcare providers through the remote health control system.</p>
57   |   | </div>
58   |   | <div class="box">
59   |   |   | 
60   |   |   | <div class="content">
61   |   |   |   | 
62   |   |   |   | <h3>Medication reminders and management</h3>
63   |   |   |   | <p>The remote health control system can provide medication reminders to help patients adhere to their medication schedules. Additionally, it can remind patients about upcoming appointments and prescription renewals.</p>
64   |   | </div>
65   |   | <div class="box">
66   |   |   | 
67   |   |   | <div class="content">
68   |   |   |   | 
69   |   |   |   | <h3>Health education and guidance</h3>
70   |   |   |   | <p>The remote health control platform can offer relevant health education and guidance to patients, including disease management, lifestyle advice, and self-care tips. It can provide access to health articles, videos, and interactive tools to support patient empowerment and informed decision-making.</p>
71   |   | </div>
72   |   | <div class="box">
73   |   |   | 
74   |   |   | <div class="content">
75   |   |   |   | 
76   |   |   |   | <h3>Emergency calls and rescue</h3>
77   |   |   |   | <p>In emergency situations, patients can use the remote health control system to send emergency calls or distress signals to promptly receive medical attention and support from emergency services.</p>
78   |   | </div>
79   |   | </div>
80   |   | </div>
81   | </div>
82   | </div>
83   | </div>
84 </section>
85
86 <section class="benefit">
87   | <div class="box">
88   |   | 
89   |   | <h3>Create Your Own Profile</h3>
90   | </div>
91   | <div class="box">
92   |   | 
93   |   | <h3>Upload Your Body Values</h3>
94   | </div>
95   | <div class="box">
96   |   | 
97   |   | <h3>Consult A Doctor Online</h3>
98   | </div>
99
100  |   | <div class="box">
101  |   |   | 
102  |   |   | <h3>Reserve</h3>
103  |   | </div>
104  | </section>
105  | <script src="Home.js"></script>
106 </body>
107 </html>
108

```

Home.js

```
① 1 let menu = document.querySelector('#Services-bar');
 2 let navbar = document.querySelector('.navbar');
 3
 4 ▼ menu.onclick = () =>{
 5   menu.classList.toggle('fa-times');
 6   navbar.classList.toggle('active');
 7 }
 8
 9 ▼ window.onscroll = () =>{
10   menu.classList.remove('fa-times');
11   navbar.classList.remove('active');
12 ▼ if(window.scrollY > 60){
13   | document.querySelector('#scroll-top').classList.add('active');
14 ▼ }else{
15   | document.querySelector('#scroll-top').classList.remove('active');
16 }
17 }
18
19 // 获取搜索按钮和输入框的引用
20 var searchButton = document.getElementById("searchButton");
21 var searchInput = document.getElementById("searchInput");
22
23 // 添加点击事件监听器
24 ▼ searchButton.addEventListener("click", function(event) {
25   event.preventDefault(); // 阻止表单提交
26   var searchText = searchInput.value;
27   console.log("搜索内容：" + searchText);
28   // 可以在这里执行搜索相关的操作
29 });
30
31 const imageLinks = document.querySelectorAll('.image-link');
32
33 ▼ imageLinks.forEach(function(link) {
34 ▼ link.addEventListener('click', function(event) {
35   | event.preventDefault();
36   | window.location.href = this.getAttribute('href');
37   });
38 });


```

Home.css

```
1 @import url('https://fonts.googleapis.com/css2?family=Agbalumo&display=swap');
```

```
2
```

```
3 ▼ :root{
```

```
4   --green:#306754;
```

```
5 }
```

```
6
```

```
7 ▼ *{
```

```
8   font-family: 'Agbalumo', Regular 400;
```

```
9   margin: 0;
```

```
10  padding: 0;
```

```
11  box-sizing: border-box;
```

```
12  outline: none;
```

```
13  border: none;
```

```
14  text-decoration: none;
```

```
15  text-transform: capitalize;
```

```
16  transition: all .2s linear;
```

```
17 }
```

```
18
```

```
19 ▼ ::selection{
```

```
20   background:var(--green);
```

```
21   color:#fff;
```

```
22 }
```

```
23 ▼ section{
```

```
24   padding:2rem 9%;
```

```
25 }
```

```
26
```

```
27 ▼ html{
```

```
28   font-size: 62.5%;
```

```
29   overflow-x: hidden;
```

```
30   scroll-behavior: smooth;
```

```
31   scroll-padding-top: 6rem;
```

```
32 }
```

```
33
```

```
34 ▼ body{
```

```
35   background:#f7f7f7;
```

```
36 }
```

```
37
```

```
38 ▼ header{
```

```
39   position: fixed;
```

```
40   top:0; left: 0; right:0;
```

```
41   z-index: 1000;
```

```
42   display: flex;
```

```
43   align-items: center;
```

```
44   justify-content: space-between;
```

```
45   background:#fff;
```

```
46   padding:2rem 9%;
```

```
47   box-shadow: 0 .5rem 1rem rgba(0,0,0,.1);
```

```
48 }
```

```
49
```

```
50 ▼ header .logo{
```

```
51   font-size: 2.5rem;
```

```
52   font-weight: bolder;
```

```
53   color:#666:
```

```
54  }
55
56 ▼ header .navbar a{
57   font-size: 2rem;
58   margin-left: 2rem;
59   color:#666;
60 }
61 ▼ header .navbar a:hover{
62   color:lightblue;
63 }
64
65 /*home*/
66 ▼ #menu-bar{
67   font-size: 3rem;
68   cursor: pointer;
69   color:#666;
70   border:.1rem solid #666;
71   border-radius: .3rem;
72   padding:.5rem 1.5rem;
73   display: none;
74 }
75
76 ▼ .btn{
77   display: inline-block;
78   padding:.8rem 3rem;
79   border:.2rem solid var(--green);
80   color:var(--green);
81   cursor: pointer;
82   font-size: 1.7rem;
83   border-radius: .5rem;
84   position: relative;
85   overflow: hidden;
86   z-index: 0;
87   margin-top: 1rem;
88 }
89 ▼ .btn::before{
90   content: '';
91   position: absolute;
92   top:0; right: 0;
93   width:0%;
94   height:100%;
95   background:var(--green);
96   transition: .3s linear;
97   z-index: -1;
98 }
99 ▼ .btn:hover::before{
100  width:100%;
101  left: 0;
102 }
103 ▼ .btn:hover{
104  color:#fff;
105 }
106
```

```
107 ▼ .home{
108   |   display: flex;
109  |   flex-wrap: wrap;
110  |   grid-gap:1.5rem;
111  |   min-height: 100vh;
112  |   align-items: center;
113  |   background:url(..../images/home-bg.jpg) no-repeat;
114  |   background-size: cover;
115  |   background-position: center;
116  |
117 ▼ .home .content{
118   |   flex:1 1 40rem;
119   |   padding-top: 6.5rem;
120  |
121 ▼ .home .image{
122   |   flex:1 1 40rem;
123  }
124 ▼ .home .image img{
125   |   width:100%;
126   |   padding:1rem;
127  }
128
129 ▼ @keyframes float{
130   |   0%, 100%{
131   |   |   transform: translateY(0rem);
132   |   }
133   |   50%{
134   |   |   transform: translateY(3rem);
135   |   }
136  }
137
138 ▼ .home .content h3{
139   |   font-size: 5rem;
140   |   color:#333;
141  }
142 ▼ .home .content p{
143   |   font-size: 1.7rem;
144   |   color:#666;
145   |   padding:1rem 0;
146  }
147 /*menu*/
148 ▼ .heading{
149   |   text-align: center;
150   |   font-size: 3.5rem;
151   |   padding:1rem;
152   |   color:#666;
153  }
154 ▼ .heading span{
155   |   color:var(--green);
156  }
157
158 ▼ .Services .box-container{
159   |   display: flex;
```

```

160 |   flex-wrap: wrap;
161 |   grid-gap: 1.5rem;
162 |
163 ▼ .Serviecs .box-container .box{
164 |   flex:1 1 30rem;
165 |   position: relative;
166 |   overflow: hidden;
167 |   box-shadow: 0 .5rem 1rem rgba(0,0,0,.1);
168 |   border:.1rem solid rgba(0,0,0,.3);
169 |   cursor: pointer;
170 |   border-radius: .5rem;
171 }
172 ▼ .Serviecs .box-container .box .image{
173 |   height:100%;
174 |   width:100%;
175 |   object-fit: cover;
176 |   position: absolute;
177 |   top:-100%; left:0;
178 }
179 ▼ .Serviecs .box-container .box .content{
180 |   text-align: center;
181 |   background:#fff;
182 |   padding:2rem;
183 }
184 ▼ .Serviecs .box-container .box .content img{
185 |   margin:1.5rem 0;
186 }
187 ▼ .Serviecs .box-container .box .content h3{
188 |   font-size: 2.5rem;
189 |   color:#333;
190 }
191 ▼ .Serviecs .box-container .box .content p{
192 |   font-size: 1.6rem;
193 |   color:#666;
194 |   padding:1rem 0;
195 }
196 ▼ .Serviecs .box-container .box:hover .image{
197 |   top:0;
198 }
199 ▼ .Serviecs .box-container .box:hover .content{
200 |   transform: translateY(100%);
201 }
202
203 /*benefit*/
204 ▼ .benifit{
205 |   display: flex;
206 |   flex-wrap: wrap;
207 |   gap:1.5rem;
208 |   padding:1rem;
209 }
210 ▼ .benifit .box{
211 |   flex:1 1 25rem;
212 |   padding:1rem;
213 |   padding:1rem;
214 |   text-align: center;
215 }
216 ▼ .benifit .box img{
217 |   border-radius: 50%;
218 |   border:.5rem solid #fff;
219 |   box-shadow: 0.5rem 1rem rgba(0,0,0,.1);
220 }
221 ▼ .benifit .box h3{
222 |   font-size: 3rem;
223 |   color:#333;
224 |   padding:1rem;

```

Login

Customer.php

```
1 <!DOCTYPE html>
2 <html lang="en">
3
4 <head>
5   <meta charset="UTF-8">
6   <meta name="viewport" content="width=device-width, initial-scale=1.0">
7   <link rel="stylesheet" href="https://cdnjs.cloudflare.com/ajax/libs/font-awesome/6.4.2/css/all.min.css">
8   <link rel="stylesheet" href="Customer.css">
9   <title>Customer Login Page </title>
10 <style>
11   input[type='password'] {
12     font-size: 20px; /* Increase text size for password fields */
13   }
14 </style>
15 </head>
16
17 <body>
18 <?php
19 session_start();
20 require_once('../conn.php');
21 if ($_SERVER["REQUEST_METHOD"] == "POST" && isset($_POST['signUp'])) {
22   $userName = $_POST['name'];
23   $Email = $_POST['email'];
24   $tel = $_POST['tel'];
25   $password = $_POST['pw'];
26
27   // Get the maximum userID
28   $maxIDQuery = "SELECT MAX(userID) FROM `users`";
29   $result = mysqli_query($conn, $maxIDQuery);
30   $row = mysqli_fetch_array($result);
31   $maxID = $row[0];
32
33   // Determine the new userID
34   $newUserID = ($maxID === NULL) ? 1 : $maxID + 1;
35
36   // Insert new user into the database
37   $sql = "INSERT INTO users (userID, userName, Email, tel, password) VALUES (?, ?, ?, ?, ?)";
38
39   // Prepare statement
40   $stmt = mysqli_prepare($conn, $sql);
41   // Bind parameters
42   mysqli_stmt_bind_param($stmt, 'issss', $newUserID, $userName, $Email, $tel, $password);
43   // Execute query
44   mysqli_stmt_execute($stmt) or die(mysqli_error($conn));
45
46   // Close statement
47   mysqli_stmt_close($stmt);
48   // Redirect or handle the successful registration as needed
49 }
50 ?>
51 <?php
52 require_once('../conn.php'); // Make sure this path is correct
53
```

```

54 // Check if we have a POST request
55 ▼ if ($_SERVER['REQUEST_METHOD'] == "POST" && isset($_POST['signIn'])) {
56     $userName = $_POST['name'];
57     $password = $_POST['pw'];
58
59     // Prepare SQL to prevent SQL injection
60     $sql = "SELECT * FROM users WHERE userName = ? AND password = ?";
61
62     // Prepare the statement
63     $stmt = mysqli_prepare($conn, $sql);
64
65     // Bind parameters
66     mysqli_stmt_bind_param($stmt, "ss", $userName, $password);
67
68     // Execute the statement
69     mysqli_stmt_execute($stmt);
70
71     // Store the result so we can check if the account exists in the database.
72     mysqli_stmt_store_result($stmt);
73
74 ▼ if (mysqli_stmt_num_rows($stmt) == 1) {
75     // Account exists, now we verify the password.
76     // Assuming you want to start a new session with user details
77     $_SESSION['userName'] = $userName;
78     // Redirect to user home page
79     header("Location: ../User/home.php");
80     exit();
81 ▼ } else {
82     // Username or password is incorrect
83     echo '<script>alert("Incorrect username or password!");</script>';
84 }
85
86 // Close statement
87 mysqli_stmt_close($stmt);
88 }
89 ?>
90 <h1>Remote Health Monitoring System</h1>
91 <h1>Login Page</h1>
92 ▼ <div class="container" id="container">
93 ▼   <div class="form-container sign-up">
94 ▼     <form method="POST">
95       <h1>Create Accounts</h1>
96       <span>use your username and email for registration</span>
97       <input type="text" id="name" name="name" placeholder="Username" required>
98       <input type="email" id="email" name="email" placeholder="Email" required>
99       <input type="text" id="tel" name="tel" placeholder="Tel" required>
100      <input type="password" id="pw" name="pw" placeholder="Password" pattern="^(?=.*\d)(?=.*[a-z])(?=.*[A-Z]).{8,}" title="Must contain at least one number and one uppercase and lowercase letter, and at least 8 or more characters" required>
101      <button id="sign-up" name="signUp">Sign Up</button>
102    </form>
103  </div>
104  <div class="form-container sign-in">
105    <form method="post" action="Customer.php">
106
107    <h1>Public</h1>
108    <h1>Sign In</h1>
109    <span>use your username password</span>
110    <input type="text" name="name" placeholder="Username">
111    <input type="password" name="pw" placeholder="Password">
112    <a href="#">Forgot Your Password?</a>
113    <button id="sign-in" name="signIn">Sign In</button>
114    <button>reset</button>
115  </form>
116 </div>
117 ▼ <div class="toggle-container">
118   <div class="toggle">
119     <div class="toggle-panel toggle-left">
120       <h1>Welcome Back!</h1>
121       <p>Enter your personal details</p>
122       <button class="hidden" id="login">Sign In</button>
123     </div>
124     <div class="toggle-panel toggle-right">
125       <h1>Hello, Friend!</h1>
126       <p>Register with your personal details</p>
127       <button class="hidden" id="register">Sign Up</button>
128       <a href="Doctor.php"><button class="restaurant">Doctor Login</button></a>
129     </div>
130   </div>
131 </div>
132 </div>
133 <script src="https://code.jquery.com/jquery-3.6.0.min.js"></script>
134 <script src="Customer.js"></script>
135 </body>
136 </html>

```

Doctor.php

```
1  <!DOCTYPE html>
2  <html lang="en">
3
4  <head>
5      <meta charset="UTF-8">
6      <meta name="viewport" content="width=device-width, initial-scale=1.0">
7      <link rel="stylesheet" href="https://cdnjs.cloudflare.com/ajax/libs/font-awesome/6.4.2/css/all.min.css">
8      <link rel="stylesheet" href="Customer.css">
9      <title>Customer Login Page </title>
10 </head>
11 <style>
12     input[type='password'] {
13         font-size: 20px; /* Increase text size for password fields */
14     }
15 </style>
16 </body>
17 <?php
18 session_start();
19 require_once('../conn.php');
20
21 if ($_SERVER["REQUEST_METHOD"] == "POST" && isset($_POST['signUp'])) {
22     // Collect and sanitize the input data
23     $name = mysqli_real_escape_string($conn, $_POST['name']);
24     $email = mysqli_real_escape_string($conn, $_POST['email']);
25     $tel = mysqli_real_escape_string($conn, $_POST['tel']);
26     $types = mysqli_real_escape_string($conn, $_POST['types']);
27     $password = mysqli_real_escape_string($conn, $_POST['pw']); // 您要求不使用密码哈希
28
29     // Retrieve the maximum doctorID and add 1 to it to create a new ID
30     $result = mysqli_query($conn, "SELECT MAX(doctorID) FROM doctor");
31     $row = mysqli_fetch_array($result);
32     $newId = str_pad((int) $row[0] + 1, 10, "0", STR_PAD_LEFT); // 假设doctorID是一个数字, 需要转换成字符串
33
34     // SQL query to insert data into the doctor table with the new doctorID
35     $sql = "INSERT INTO doctor (doctorID, name, types, tel, email, password) VALUES (?, ?, ?, ?, ?, ?)";
36
37     // Prepare the SQL statement
38     if ($stmt = mysqli_prepare($conn, $sql)) {
39         // Bind variables to the prepared statement as parameters
40         mysqli_stmt_bind_param($stmt, "ssssss", $newId, $name, $types, $tel, $email, $password);
41         // Execute the prepared statement
42         if (mysqli_stmt_execute($stmt)) {
43             echo "New record created successfully with ID " . $newId;
44         } else {
45             echo "Error: " . mysqli_error($conn);
46         }
47     }
48
49     // Close statement
50     mysqli_stmt_close($stmt);
51 } else {
52     echo "Error: " . mysqli_error($conn);
53 }
```

```

54      |
55      | // Close connection
56      | mysqli_close($conn);
57  }
58 ?>
59 <?php
60 require_once('../conn.php'); // Make sure this path is correct
61
62 // Check if we have a POST request
63 if ($_SERVER["REQUEST_METHOD"] == "POST" && isset($_POST['signIn'])) {
64     $userName = $_POST['name'];
65     $password = $_POST['pw'];
66
67     // Prepare SQL to prevent SQL injection
68     $sql = "SELECT * FROM doctor WHERE name = ? AND password = ?";
69
70     // Prepare the statement
71     $stmt = mysqli_prepare($conn, $sql);
72
73     // Bind parameters
74     mysqli_stmt_bind_param($stmt, "ss", $userName, $password);
75
76     // Execute the statement
77     mysqli_stmt_execute($stmt);
78
79     // Store the result so we can check if the account exists in the database.
80     mysqli_stmt_store_result($stmt);
81
82 if (mysqli_stmt_num_rows($stmt) == 1) {
83     // Account exists, now we verify the password.
84     // Assuming you want to start a new session with user details
85     $_SESSION['name'] = $userName;
86     // Redirect to user home page
87     header("Location: ../Doctor/home.php");
88     exit();
89 } else {
90     // Username or password is incorrect
91     echo '<script>alert("Incorrect username or password!");</script>';
92 }
93
94 // Close statement
95 mysqli_stmt_close($stmt);
96 }
97 ?>
98
99 <h1>Remote Health Monitoring System</h1>
100 <h1>Login Page</h1>
101 <div class="container" id="container">
102 <div class="form-container sign-up">
103 <form method="POST">
104     <h1>Create Account</h1>
105     <span>use your username and email for registration</span>
106     <input type="text" id="name" name="name" placeholder="Username" required>
107

```

```

107 |     <input type="email" id="email" name="email" placeholder="Email" required>
108 |     <input type="tel" id="tel" name="tel" placeholder="Tel" required>
109 |     <input list="types" id="types-input" name="types" placeholder="Doctor type" required/>
110 ▼    <datalist id="types">
111 |        <option value="General medical practitioner">
112 |        <option value="Cardiologist">
113 |        <option value="Gynocologist">
114 |        <option value="Pediatrician">
115 |    </datalist>
116 |    <input type="password" name="pw" placeholder="Password" pattern="(?=.+\d)(?=.[a-z])(?=.*[A-Z]).{8,}" title="Must contain at least one number and one uppercase and lowercase letter, and at least 8 or more characters" required>
117 |    <button id="sign-up" name="signUp">Sign Up</button>
118 |
119 |</form>
120 |</div>
121 ▼ <div class="form-container sign-in">
122 ▼   <form method="post" action="Doctor.php">
123 |       <h1>Doctor</h1>
124 |       <h1>Sign In</h1>
125 |       <span>use your username password</span>
126 |       <input type="text" name="name" placeholder="Username">
127 |       <input type="password" name="pw" placeholder="Password">
128 |       <a href="#">Forget Your Password?</a>
129 |       <button id="sign-in" name="signin">Sign In</button>
130 |       <button>reset</button>
131 |   </form>
132 |</div>
133 ▼ <div class="toggle-container">
134 |   <div class="toggle">
135 ▼     <div class="toggle-panel toggle-left">
136 |         <h1>Welcome Back!</h1>
137 |         <p>Enter your personal details:</p>
138 |         <button class="hidden" id="login">Sign In</button>
139 |     </div>
140 ▼     <div class="toggle-panel toggle-right">
141 |         <h1>Hello, Friend!</h1>
142 |         <p>Register with your personal details</p>
143 |         <button class="hidden" id="register">Sign Up</button>
144 |         <a href="Customer.php"><button class="restaurant">Public Login</button></a>
145 |     </div>
146 |   </div>
147 |</div>
148 |<script src="https://code.jquery.com/jquery-3.6.0.min.js"></script>
149 |<script src="Customer.js"></script>
150 |
151 |</body>
152 |</html>

```

RegisterDoctor.php

```
1  <!DOCTYPE html>
2  <html lang="en">
3
4  <head>
5      <meta charset="UTF-8">
6      <title>Register</title>
7
8      <script src="https://code.jquery.com/jquery-3.6.0.min.js"></script>
9      <script src="https://kit.fontawesome.com/61b76a306e.js" crossorigin="anonymous"></script>
10     <script src="DoctorNext.js"></script>
11     <script src="password.js"></script>
12     <script src="inputValidation.js"></script>
13     <link rel="stylesheet" href="index.css">
14     <link rel="stylesheet" href="styles.css">
15     <link rel="stylesheet" href="termsConditions.css">
16     <script src="function.js"></script>
17 </head>
18 <script>
19     function navigate() {
20         var dropdown = document.getElementById("role");
21         var selectedOption = dropdown.options[dropdown.selectedIndex].value;
22
23         if (selectedOption !== "") {
24             window.location.href = selectedOption;
25         }
26     }
27
28     // Function to toggle between registration form and password form
29     function toggleForms() {
30         var registrationForm = document.getElementById('registrationForm');
31         var passwordForm = document.getElementById('passwordForm');
32         var nextButton = document.getElementById('next');
33
34         // When the Next button is clicked
35         document.getElementById('nextButton').addEventListener('click', function(event){
36             event.preventDefault(); // Prevent the form from submitting
37
38             // Hide the registration form and show the password form
39             registrationForm.style.display = 'none';
40             passwordForm.classList.remove('hidden');
41         });
42     }
43
44     // Call the function to activate form toggling
45     toggleForms();
46
47 </script>
48 <body>
49 <?
50 require_once("../conn.php");
51 // Check if the form is submitted
52 if ($_SERVER["REQUEST_METHOD"] == "POST") {
53     // Get the maximum doctorID from the table
```

```

56     $query = "SELECT MAX(doctorID) as maxID FROM doctor";
57     $result = mysqli_query($conn, $query);
58     $row = mysqli_fetch_assoc($result);
59     $maxId = $row['maxID'];
60     $newId = ($int)$maxId + 1; // Increment the doctorID
61
62     // Get the form data
63     $name = mysqli_real_escape_string($conn, $_POST['name']);
64     $types = mysqli_real_escape_string($conn, $_POST['types']);
65     $tel = mysqli_real_escape_string($conn, $_POST['tel']);
66     $email = mysqli_real_escape_string($conn, $_POST['email']);
67     $password = mysqli_real_escape_string($conn, $_POST['password']); // Consider hashing the password
68
69     // Prepare an insert statement
70     $query = "INSERT INTO doctor (doctorID, name, types, tel, email, password) VALUES ($newId, $name, $types, $tel, $email, $password)";
71
72     if($stmt = mysqli_prepare($conn, $query)){
73         // Bind parameters to the prepared statement as parameters
74         mysqli_stmt_bind_param($stmt, "isiss", $newId, $name, $types, $tel, $email, $password);
75
76         // Execute the statement
77         if(mysqli_stmt_execute($stmt)){
78             echo "Record created successfully.";
79         } else {
80             echo "ERROR: Could not execute query: $query. " . mysqli_error($conn);
81         }
82     }
83
84     // Close statement
85     mysqli_stmt_close($stmt);
86
87     // Close connection
88     mysqli_close($conn);
89 }
90 </div>
91 <div class="form">
92 | <div>Doctor Register</div>
93 |
94 | <!-- personal information form -->
95 | <form id="registrationForm" method="post" action="RegisterDoctor.php">
96 |   <div class="field username-field">
97 |     <div class="input">
98 |       <input type="text" name="name" id="username" pattern="^([A-Z][a-z]* required placeholder="start with an uppercase letter and contain lowercase letters">
99 |     <label for="username">Username:</label>
100 |   </div>
101 |   <span class="error username-error"><i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> start from uppercase and contain lowercase</span>
102 | </div>
103 | </div>
104 | <div class="field types-field" id="types-field">
105 |   <label for="types">Doctor Type:</label>

```

```

107 |   <div class="input" name="types">
108 |     <input list="types" />
109 |     <datalist id="types">
110 |       <option value="General medical practitioner" />
111 |       <option value="Cardiologist" />
112 |       <option value="Gynecologist" />
113 |       <option value="Pediatrician" />
114 |     </datalist>
115 |
116 |   </div>
117 |
118 |   <span>
119 |     <div class="field email-password">
120 |       <div class="input">
121 |         <input type="email" name="email" id="email" required>
122 |         <label for="email">Email:</label>
123 |       </div>
124 |       <span class="error email-error">
125 |         <p class="error-text"> <i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> Please enter a valid email</p>
126 |       </span>
127 |     </div>
128 |
129 |     <div class="field tel-field">
130 |       <div class="input">
131 |         <input type="text" name="tel" id="tel" required pattern="^852-\d{8}" value="852->
132 |         <label for="tel">Tel:</label>
133 |       </div>
134 |       <span class="error tel-error">
135 |         <p class="error-text"> <i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> tel format: 852- 8 numbers.</p>
136 |       </span>
137 |     </div>
138 |
139 |     <div class="field role-field">
140 |       <div class="input">
141 |         User Type:
142 |         <select id="role" onchange="navigate()">
143 |           <option value="" disabled selected>Select your role</option>
144 |           <option value="RegisterUser.php">Public</option>
145 |           <option value="Doctor" selected>Doctor</option>
146 |         </select>
147 |       </div>
148 |       <span class="error role-error">
149 |         <p class="error-text"> <i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> Please select a role</p>
150 |       </span>
151 |     </div>
152 |     already have an account? <a href="../index.html">Login</a>
153 |
154 |     <button type="button" id="next"><i class="fa-solid fa-arrow-right"></i> Next</button>
155 |   </form>
156 |
157 | <!-- password form -->
158 | <form id="passwordForm" class="hidden" method="post" action="RegisterDoctor.php">

```

```

160      <div class="field create-password">
161        <div class="input">
162          <input type="password" name="password" id="password" minlength="8" required>
163          <label for="password">Password:</label>
164          <span class="eye" onclick="Actpassword('password')">
165            <i id="password-hide" class="fa-regular fa-eye-slash"></i>
166            <i id="password-show" class="fa-regular fa-eye" style="display: none;"></i>
167          </span>
168        </div>
169
170        <span class="error password-error">
171          <i class="bx bx-error-circle error-icon"></i>
172          <p class="Starterror-text">
173            <i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> Start with an uppercase letter.
174          </p>
175          <p class="Minerror-text">
176            <i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> Minimum one lowercase letter and one digit.
177          </p>
178          <p class="Lengtherror-text">
179            <i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> Length must be between 8 to 10 characters.
180          </p>
181        </span>
182      </div>
183
184      <div class="field confirm-password">
185        <div class="input">
186          <input type="password" id="confirmPassword" required>
187          <label for="confirmPassword">Confirm Password:</label>
188          <span class="eye" onclick="Actpassword('confirmPassword')">
189            <i id="confirmPassword-hide" class="fa-regular fa-eye-slash"></i>
190            <i id="confirmPassword-show" class="fa-regular fa-eye" style="display: none;"></i>
191          </span>
192        </div>
193
194        <div id="passwordError" class="error-message">Passwords do not match.</div>
195        <span class="cPassword-error">
196          <p class="error-text"> <i class="fa-solid fa-circle-exclamation" style="color: #e61405;"></i> Passwords don't match</p>
197        </span>
198      </div>
199      <input type="checkbox" required>I agree to the <a id="termsLink" href="#">terms and conditions</a>
200      <br>
201      <button id="back"><i class="fa-solid fa-arrow-left"></i> Back</button>
202      <button type="submit" id="register">Register</button>
203    </form>
204  </div>
205
206  <!-- Terms and conditions overlay -->
207  <div class="overlay">
208    <div class="modal">
209      <span class="modal-close">&times;;</span>
210      <h2>Terms and Conditions</h2>
211      <p>1. Account Registrations</p>
212    </div>

```

```

213      <ul style="list-style-type: none; padding-left: 0; margin: 0; border-left: 1px solid black; padding-left: 10px; margin-left: 20px;">
214        <li>1.1. To use our food delivery platform, you must register an account with accurate and complete information.</li>
215        <li>1.2. You are responsible for maintaining the confidentiality of your account credentials and for all activities that occur under your account.<>
216        <li>1.3. You must promptly notify us of any unauthorized use of your account or any other security breach.</li>
217        <li>2. Placing Orders</li>
218        <li>2.1. By placing an order through our platform, you agree to provide accurate and up-to-date information, including your delivery address and payment method.
219        <li>2.2. You are responsible for reviewing and confirming your order details before submitting the order.</li>
220        <li>2.3. Once an order is placed, it is considered final and cannot be modified or canceled, except as permitted by our cancellation policy.
221
222        <li>3. Delivery</li>
223        <li>3.1. We strive to deliver orders promptly and within the estimated delivery time. However, delivery times may vary depending on factors such as traffic, weather, and order volume.</li>
224        <li>3.2. You agree to provide clear and accurate delivery instructions to ensure smooth delivery of your order.</li>
225        <li>3.3. In the event of any issues or delays in delivery, we will make reasonable efforts to inform you and resolve the situation as quickly as possible.
226
227        <li>4. Payment</li>
228        <li>4.1. Payment for orders placed through our platform must be made using the available payment methods.</li>
229        <li>4.2. All prices displayed on our platform are in the local currency and include applicable taxes and fees, unless otherwise specified.
230        <li>4.3. We reserve the right to update prices, fees, and payment methods at any time without prior notice.</li>
231
232        <li>5. User Conduct</li>
233        <li>5.1. You agree to use our food delivery platform responsibly and in compliance with all applicable laws and regulations.
234        <li>5.2. You shall not engage in any fraudulent, unlawful, or abusive activities while using our services.</li>
235
236        <li>5.3. Any misuse or unauthorized use of our platform may result in the suspension or termination of your account.</li>
237        <li>6. Intellectual Property</li>
238        <li>6.1. All content, trademarks, logos, and intellectual property rights on our platform are owned or licensed by us.</li>
239        <li>6.2. You may not reproduce, distribute, modify, or create derivative works of any content from our platform without our permission.
240
241        <li>7. Limitation of Liability</li>
242        <li>7.1. We strive to provide accurate information and a seamless experience. However, we do not guarantee the availability, reliability, or performance of our services.
243        <li>7.2. In no event shall we be liable for any direct, indirect, incidental, consequential, or punitive damages arising out of the use of our platform.
244
245        <li>8. Modifications to Terms and Conditions
246        <li>8.1. We reserve the right to modify these terms and conditions at any time without prior notice.</li>
247        <li>8.2. It is your responsibility to review these terms periodically for any changes.</li>
248        <li>By using our food delivery platform, you acknowledge that you have read, understood, and agreed to these terms and conditions. If you have any questions or concerns, please contact our customer support.</li>
249
250        <li>Please note that this is just a sample terms and conditions content, and it's important to consult with legal professionals to ensure it complies with applicable laws and regulations.</li>
251
252    </ul>
253  </div>
254
255</body>
256</html>
257
258
259

```

Customer.js

```
 1  const container = document.getElementById('container');
 2  const registerBtn = document.getElementById('register');
 3  const loginBtn = document.getElementById('login');
 4
 5 ▼ registerBtn.addEventListener('click', () => {
 6  |   container.classList.add("active");
 7  });
 8
 9 ▼ loginBtn.addEventListener('click', () => {
10  |   container.classList.remove("active");
11  });
12
13 ▼ $(document).ready(function() {
14
15
16  |   $('#reset').click(function() {
17  |     $('#username').val('');
18  |     $('#welcome-message').empty();
19  |   });
20
21  |   $('#sign-up').click(function() {
22  |     var name = $('#name').val();
23  |     if (name !== '') {
24  |       alert('You have successfully registered. Please verify your email.');
25  |     }
26  |   });
27  });
28
29
30 ▼ $(document).ready(function() {
31  |   $('#login').addClass('striped');
32
33  |   $('tr').hover(
34  |     function() {
35  |       $(this).addClass('over');
36  |     },
37  |     function() {
38  |       $(this).removeClass('over');
39  |     }
40  |   );
41  });
42
```

Customer.js

```
1  @import url('https://fonts.googleapis.com/css2?family=Montserrat:wght@300;400;500;600;700&display=swap');
2
3  *{
4      margin: 0;
5      padding: 0;
6      box-sizing: border-box;
7      font-family: 'Montserrat', sans-serif;
8  }
9
10 body{
11     background-color: #c9d6ff;
12     background: linear-gradient(to right, #e2e2e2, #c9d6ff);
13     display: flex;
14     align-items: center;
15     justify-content: center;
16     flex-direction: column;
17     height: 100vh;
18 }
19
20 .container{
21     background-color: #fff;
22     border-radius: 30px;
23     box-shadow: 0 5px 15px rgba(0, 0, 0, 0.35);
24     position: relative;
25     overflow: hidden;
26     width: 768px;
27     max-width: 100%;
28     min-height: 480px;
29 }
30
31 .container p{
32     font-size: 14px;
33     line-height: 20px;
34     letter-spacing: 0.3px;
35     margin: 20px 0;
36 }
37
38 .container span{
39     font-size: 12px;
40 }
41
42 .container a{
43     color: #333;
44     font-size: 13px;
45     text-decoration: none;
46     margin: 15px 0 10px;
47 }
48
49 .container button{
50     background-color: #512da8;
51     color: #fff;
52     font-size: 12px;
53     padding: 10px 45px;
```

```
54 |     border: 1px solid transparent;
55 |     border-radius: 8px;
56 |     font-weight: 600;
57 |     letter-spacing: 0.5px;
58 |     text-transform: uppercase;
59 |     margin-top: 10px;
60 |     cursor: pointer;
61 }
62
63 ▼ .container button.hidden{
64 |     background-color: transparent;
65 |     border-color: #fff;
66 }
67
68 ▼ .container button.restaurant{
69 |     background-color: rgb(3, 70, 3);
70 |     border-color: #fff;
71 }
72
73 ▼ .container button.delivery{
74 |     background-color: rgb(104, 4, 4);
75 |     border-color: #fff;
76 }
77
78 ▼ .container form{
79 |     background-color: #fff;
80 |     display: flex;
81 |     align-items: center;
82 |     justify-content: center;
83 |     flex-direction: column;
84 |     padding: 0 40px;
85 |     height: 100%;
86 }
87
88 ▼ .container input{
89 |     background-color: #eee;
90 |     border: none;
91 |     margin: 8px 0;
92 |     padding: 10px 15px;
93 |     font-size: 13px;
94 |     border-radius: 8px;
95 |     width: 100%;
96 |     outline: none;
97 }
98
99 ▼ .container button:hover {
100 |     background-color: #80bfff;
101 }
102
103 ▼ .form-container{
104 |     position: absolute;
105 |     top: 0;
106 |     height: 100%;
```

```
107 |   transition: all 0.6s ease-in-out;
108 }
109
110 ▼ .sign-in{
111   left: 0;
112   width: 50%;
113   z-index: 2;
114 }
115
116 ▼ .container.active .sign-in{
117   transform: translateX(100%);
118 }
119
120 ▼ .sign-up{
121   left: 0;
122   width: 50%;
123   opacity: 0;
124   z-index: 1;
125 }
126
127 ▼ .container.active .sign-up{
128   transform: translateX(100%);
129   opacity: 1;
130   z-index: 5;
131   animation: move 0.6s;
132 }
133
134 ▼ @keyframes move{
135   0%, 49.99%{
136     opacity: 0;
137     z-index: 1;
138   }
139   50%, 100%{
140     opacity: 1;
141     z-index: 5;
142   }
143 }
144
145 ▼ .toggle-container{
146   position: absolute;
147   top: 0;
148   left: 50%;
149   width: 50%;
150   height: 100%;
151   overflow: hidden;
152   transition: all 0.6s ease-in-out;
153   border-radius: 150px 0 0 100px;
154   z-index: 1000;
155 }
156
157 ▼ .container.active .toggle-container{
158   transform: translateX(-100%);
159   border-radius: 0 150px 100px 0;
160 }
```

```
161
162 ▼ .toggle{
163     background-color: #512da8;
164     height: 100%;
165     background: linear-gradient(to right, #5c6bc0, #512da8);
166     color: #fff;
167     position: relative;
168     left: -100%;
169     height: 100%;
170     width: 200%;
171     transform: translateX(0);
172     transition: all 0.6s ease-in-out;
173 }
174
175 ▼ .container.active .toggle{
176     transform: translateX(50%);
177 }
178
179 ▼ .toggle-panel{
180     position: absolute;
181     width: 50%;
182     height: 100%;
183     display: flex;
184     align-items: center;
185     justify-content: center;
186     flex-direction: column;
187     padding: 0 30px;
188     text-align: center;
189     top: 0;
190     transform: translateX(0);
191     transition: all 0.6s ease-in-out;
192 }
193
194 ▼ .toggle-left{
195     transform: translateX(-200%);
196 }
197
198 ▼ .container.active .toggle-left{
199     transform: translateX(0);
200 }
201
202 ▼ .toggle-right{
203     right: 0;
204     transform: translateX(0);
205 }
206
207 ▼ .container.active .toggle-right{
208     transform: translateX(200%);
209 }
```

Meeting Room

Meeting Room.html

```
1  <!DOCTYPE html>
2  <html lang="en">
3
4  <head>
5      <meta charset="UTF-8">
6      <title>Document</title>
7
8  </head>
9
10 <body>
11 <script>
12     var script = document.createElement("script");
13     script.type = "text/javascript";
14
15     script.addEventListener("load", function(event) {
16         const config = {
17             name: "",
18             meetingId: "fyp",
19             apiKey: "34d073b0-96b3-44df-8a26-7d2450faf7e3",
20
21             containerId: null,
22
23             micEnabled: true,
24             webcamEnabled: true,
25             participantCanToggleSelfWebcam: true,
26             participantCanToggleSelfMic: true,
27
28             chatEnabled: true,
29             raiseHandEnabled: true,
30             screenShareEnabled: true,
31             whiteboardEnabled: true,
32
33             joinScreen: {
34                 visible: true,
35             },
36
37             recording: {
38                 enabled: true,
39                 webhookUrl: "https://www.videosdk.live/callback",
40                 autoStart: false,
41                 theme: "DARK", // DARK || LIGHT || DEFAULT
42
43             layout: {
44                 type: "SIDEBAR",
45                 priority: "PIN",
46                 gridSize: 3,
47             },
48         },
49
50         permissions: {
51             pin: true,
52             toggleRecording: true,
53             askToJoin: false,
54             toggleParticipantWebcam: false,
55             toggleParticipantMic: false,
56             toggleParticipantScreenshare: false,
57             removeParticipant: true,
58             endMeeting: true,
59             drawOnWhiteboard: true,
60             toggleWhiteboard: true,
61         },
62     });
63
64
65
66
67
68
69
70
71
72
73
74
75
76
77
78
79
80
81
82
83
84
85
86
87
88
89
90
91
92
93
94
95
96
97
98
99
100
101
102
103
104
105
106
107
108
109
110
111
112
113
114
115
116
117
118
119
120
121
122
123
124
125
126
127
128
129
130
131
132
133
134
135
136
137
138
139
140
141
142
143
144
145
146
147
148
149
150
151
152
153
154
155
156
157
158
159
160
161
162
163
164
165
166
167
168
169
170
171
172
173
174
175
176
177
178
179
180
181
182
183
184
185
186
187
188
189
190
191
192
193
194
195
196
197
198
199
200
201
202
203
204
205
206
207
208
209
210
211
212
213
214
215
216
217
218
219
220
221
222
223
224
225
226
227
228
229
230
231
232
233
234
235
236
237
238
239
240
241
242
243
244
245
246
247
248
249
250
251
252
253
254
255
256
257
258
259
259
260
261
262
263
264
265
266
267
268
269
270
271
272
273
274
275
276
277
278
279
279
280
281
282
283
284
285
286
287
288
289
289
290
291
292
293
294
295
296
297
298
299
299
300
301
302
303
304
305
306
307
308
309
309
310
311
312
313
314
315
316
317
318
319
319
320
321
322
323
324
325
326
327
328
329
329
330
331
332
333
334
335
336
337
338
339
339
340
341
342
343
344
345
346
347
348
349
349
350
351
352
353
354
355
356
357
358
359
359
360
361
362
363
364
365
366
367
368
369
369
370
371
372
373
374
375
376
377
378
379
379
380
381
382
383
384
385
386
387
388
389
389
390
391
392
393
394
395
396
397
398
399
399
400
401
402
403
404
405
406
407
408
409
409
410
411
412
413
414
415
416
417
418
419
419
420
421
422
423
424
425
426
427
428
429
429
430
431
432
433
434
435
436
437
438
439
439
440
441
442
443
444
445
446
447
448
449
449
450
451
452
453
454
455
456
457
458
459
459
460
461
462
463
464
465
466
467
468
469
469
470
471
472
473
474
475
476
477
478
479
479
480
481
482
483
484
485
486
487
488
489
489
490
491
492
493
494
495
496
497
498
499
499
500
501
502
503
504
505
506
507
508
509
509
510
511
512
513
514
515
516
517
518
519
519
520
521
522
523
524
525
526
527
528
529
529
530
531
532
533
534
535
536
537
538
539
539
540
541
542
543
544
545
546
547
548
549
549
550
551
552
553
554
555
556
557
558
559
559
560
561
562
563
564
565
566
567
568
569
569
570
571
572
573
574
575
576
577
578
579
579
580
581
582
583
584
585
586
587
588
589
589
590
591
592
593
594
595
596
597
598
599
599
600
601
602
603
604
605
606
607
608
609
609
610
611
612
613
614
615
616
617
618
619
619
620
621
622
623
624
625
626
627
628
629
629
630
631
632
633
634
635
636
637
638
639
639
640
641
642
643
644
645
646
647
648
649
649
650
651
652
653
654
655
656
657
658
659
659
660
661
662
663
664
665
666
667
668
669
669
670
671
672
673
674
675
676
677
678
679
679
680
681
682
683
684
685
686
687
688
689
689
690
691
692
693
694
695
696
697
698
699
699
700
701
702
703
704
705
706
707
708
709
709
710
711
712
713
714
715
716
717
718
719
719
720
721
722
723
724
725
726
727
728
729
729
730
731
732
733
734
735
736
737
738
739
739
740
741
742
743
744
745
746
747
748
749
749
750
751
752
753
754
755
756
757
758
759
759
760
761
762
763
764
765
766
767
768
769
769
770
771
772
773
774
775
776
777
778
779
779
780
781
782
783
784
785
786
787
788
789
789
790
791
792
793
794
795
796
797
797
798
799
799
800
801
802
803
804
805
806
807
808
809
809
810
811
812
813
814
815
816
817
818
819
819
820
821
822
823
824
825
826
827
828
829
829
830
831
832
833
834
835
836
837
838
839
839
840
841
842
843
844
845
846
847
848
849
849
850
851
852
853
854
855
856
857
858
859
859
860
861
862
863
864
865
866
867
868
869
869
870
871
872
873
874
875
876
877
878
879
879
880
881
882
883
884
885
886
887
888
888
889
889
890
891
892
893
894
895
896
897
898
898
899
899
900
901
902
903
904
905
906
907
908
909
909
910
911
912
913
914
915
916
917
918
919
919
920
921
922
923
924
925
926
927
928
929
929
930
931
932
933
934
935
936
937
938
939
939
940
941
942
943
944
945
946
947
948
949
949
950
951
952
953
954
955
956
957
958
959
959
960
961
962
963
964
965
966
967
968
969
969
970
971
972
973
974
975
976
977
978
979
979
980
981
982
983
984
985
986
987
987
988
989
989
990
991
992
993
994
995
996
997
998
999
999
1000
1000
1001
1002
1003
1004
1005
1006
1007
1008
1009
1009
1010
1011
1012
1013
1014
1015
1016
1017
1018
1019
1019
1020
1021
1022
1023
1024
1025
1026
1027
1028
1029
1029
1030
1031
1032
1033
1034
1035
1036
1037
1038
1039
1039
1040
1041
1042
1043
1044
1045
1046
1047
1048
1048
1049
1050
1051
1052
1053
1054
1055
1056
1057
1058
1059
1059
1060
1061
1062
1063
1064
1065
1066
1067
1068
1069
1069
1070
1071
1072
1073
1074
1075
1076
1077
1078
1078
1079
1080
1081
1082
1083
1084
1085
1086
1087
1087
1088
1089
1089
1090
1091
1092
1093
1094
1095
1096
1096
1097
1098
1099
1099
1100
1101
1102
1103
1104
1105
1106
1107
1108
1109
1109
1110
1111
1112
1113
1114
1115
1116
1117
1118
1119
1119
1120
1121
1122
1123
1124
1125
1126
1127
1128
1129
1129
1130
1131
1132
1133
1134
1135
1136
1137
1138
1139
1139
1140
1141
1142
1143
1144
1145
1146
1147
1148
1148
1149
1150
1151
1152
1153
1154
1155
1156
1157
1158
1159
1159
1160
1161
1162
1163
1164
1165
1166
1167
1168
1169
1169
1170
1171
1172
1173
1174
1175
1176
1177
1178
1178
1179
1180
1181
1182
1183
1184
1185
1186
1187
1187
1188
1189
1189
1190
1191
1192
1193
1194
1195
1195
1196
1197
1198
1199
1199
1200
1201
1202
1203
1204
1205
1206
1207
1208
1209
1209
1210
1211
1212
1213
1214
1215
1216
1217
1218
1219
1219
1220
1221
1222
1223
1224
1225
1226
1227
1228
1229
1229
1230
1231
1232
1233
1234
1235
1236
1237
1238
1239
1239
1240
1241
1242
1243
1244
1245
1246
1247
1248
1248
1249
1250
1251
1252
1253
1254
1255
1256
1257
1258
1259
1259
1260
1261
1262
1263
1264
1265
1266
1267
1268
1269
1269
1270
1271
1272
1273
1274
1275
1276
1277
1278
1278
1279
1280
1281
1282
1283
1284
1285
1286
1287
1287
1288
1289
1289
1290
1291
1292
1293
1294
1295
1296
1297
1297
1298
1299
1299
1300
1301
1302
1303
1304
1305
1306
1307
1308
1309
1309
1310
1311
1312
1313
1314
1315
1316
1317
1318
1318
1319
1320
1321
1322
1323
1324
1325
1326
1327
1328
1329
1329
1330
1331
1332
1333
1334
1335
1336
1337
1338
1339
1339
1340
1341
1342
1343
1344
1345
1346
1347
1348
1348
1349
1350
1351
1352
1353
1354
1355
1356
1357
1358
1359
1359
1360
1361
1362
1363
1364
1365
1366
1367
1368
1369
1369
1370
1371
1372
1373
1374
1375
1376
1377
1378
1378
1379
1380
1381
1382
1383
1384
1385
1386
1387
1387
1388
1389
1389
1390
1391
1392
1393
1394
1395
1396
1396
1397
1398
1399
1399
1400
1401
1402
1403
1404
1405
1406
1407
1408
1409
1409
1410
1411
1412
1413
1414
1415
1416
1417
1418
1418
1419
1420
1421
1422
1423
1424
1425
1426
1427
1428
1428
1429
1430
1431
1432
1433
1434
1435
1436
1437
1438
1438
1439
1440
1441
1442
1443
1444
1445
1446
1447
1448
1448
1449
1450
1451
1452
1453
1454
1455
1456
1457
1458
1459
1459
1460
1461
1462
1463
1464
1465
1466
1467
1468
1468
1469
1470
1471
1472
1473
1474
1475
1476
1477
1478
1478
1479
1480
1481
1482
1483
1484
1485
1486
1487
1488
1488
1489
1490
1491
1492
1493
1494
1495
1496
1496
1497
1498
1499
1499
1500
1501
1502
1503
1504
1505
1506
1507
1508
1509
1509
1510
1511
1512
1513
1514
1515
1516
1517
1518
1518
1519
1520
1521
1522
1523
1524
1525
1526
1527
1528
1528
1529
1530
1531
1532
1533
1534
1535
1536
1537
1538
1538
1539
1540
1541
1542
1543
1544
1545
1546
1547
1548
1548
1549
1550
1551
1552
1553
1554
1555
1556
1557
1558
1559
1559
1560
1561
1562
1563
1564
1565
1566
1567
1568
1568
1569
1570
1571
1572
1573
1574
1575
1576
1577
1578
1578
1579
1580
1581
1582
1583
1584
1585
1586
1587
1588
1588
1589
1590
1591
1592
1593
1594
1595
1596
1596
1597
1598
1599
1599
1600
1601
1602
1603
1604
1605
1606
1607
1608
1609
1609
1610
1611
1612
1613
1614
1615
1616
1617
1618
1618
1619
1620
1621
1622
1623
1624
1625
1626
1627
1628
1628
1629
1630
1631
1632
1633
1634
1635
1636
1637
1638
1638
1639
1640
1641
1642
1643
1644
1645
1646
1647
1648
1648
1649
1650
1651
1652
1653
1654
1655
1656
1657
1658
1659
1659
1660
1661
1662
1663
1664
1665
1666
1667
1668
1668
1669
1670
1671
1672
1673
1674
1675
1676
1677
1678
1678
1679
1680
1681
1682
1683
1684
1685
1686
1687
1688
1688
1689
1690
1691
1692
1693
1694
1695
1696
1696
1697
1698
1699
1699
1700
1701
1702
1703
1704
1705
1706
1707
1708
1709
1709
1710
1711
1712
1713
1714
1715
1716
1717
1718
1718
1719
1720
1721
1722
1723
1724
1725
1726
1727
1728
1728
1729
1730
1731
1732
1733
1734
1735
1736
1737
1738
1738
1739
1740
1741
1742
1743
1744
1745
1746
1747
1748
1748
1749
1750
1751
1752
1753
1754
1755
1756
1757
1758
1759
1759
1760
1761
1762
1763
1764
1765
1766
1767
1768
1768
1769
1770
1771
1772
1773
1774
1775
1776
1777
1778
1778
1779
1780
1781
1782
1783
1784
1785
1786
1787
1788
1788
1789
1790
1791
1792
1793
1794
1795
1796
1796
1797
1798
1799
1799
1800
1801
1802
1803
1804
1805
1806
1807
1808
1809
1809
1810
1811
1812
1813
1814
1815
1816
1817
1818
1818
1819
1820
1821
1822
1823
1824
1825
1826
1827
1828
1828
1829
1830
1831
1832
1833
1834
1835
1836
1837
1838
1838
1839
1840
1841
1842
1843
1844
1845
1846
1847
1848
1848
1849
1850
1851
1852
1853
1854
1855
1856
1857
1858
1859
1859
1860
1861
1862
1863
1864
1865
1866
1867
1868
1868
1869
1870
1871
1872
1873
1874
1875
1876
1877
1877
1878
1879
1880
1881
1882
1883
1884
1885
1886
1887
1888
1888
1889
1890
1891
1892
1893
1894
1895
1896
1896
1897
1898
1899
1899
1900
1901
1902
1903
1904
1905
1906
1907
1908
1909
1909
1910
1911
1912
1913
1914
1915
1916
1917
1918
1918
1919
1920
1921
1922
1923
1924
1925
1926
1927
1928
1928
1929
1930
1931
1932
1933
1934
1935
1936
1937
1938
1938
1939
1940
1941
1942
1943
1944
1945
1946
1947
1948
1948
1949
1950
1951
1952
1953
1954
1955
1956
1957
1958
1959
1959
1960
1961
1962
1963
1964
1965
1966
1967
1968
1968
1969
1970
1971
1972
1973
1974
1975
1976
1977
1977
1978
1979
1980
1981
1982
1983
1984
1985
1986
1987
1988
1988
1989
1990
1991
1992
1993
1994
1995
1996
1996
1997
1998
1999
1999
2000
2001
2002
2003
2004
2005
2006
2007
2008
2009
2009
2010
2011
2012
2013
2014
2015
2016
2017
2018
2018
2019
2020
2021
2022
2023
2024
2025
2026
2027
2028
2029
2029
2030
2031
2032
2033
2034
2035
2036
2037
2038
2038
2039
2040
2041
2042
2043
2044
2045
2046
2047
2048
2048
2049
2050
2051
2052
2053
2054
2055
2056
2057
2058
2059
2059
2060
2061
2062
2063
2064
2065
2066
2067
2068
2068
2069
2070
2071
2072
2073
2074
2075
2076
2077
2078
2078
2079
2080
2081
2082
2083
2084
2085
2086
2087
2088
2088
2089
2090
2091
2092
2093
2094
2095
2096
2096
2097
2098
2099
2099
2100
2101
2102
2103
2104
2105
2106
2107
2108
2109
2109
2110
2111
2112
2113
2114
2115
2116
2117
2118
2118
2119
2120
2121
2122
2123
2124
2125
2126
2127
2128
2129
2129
2130
2131
2132
2133
2134
2135
2136
2137
2138
2138
2139
2140
2141
2142
2143
2144
2145
2146
2147
2148
2148
2149
2150
2151
2152
2153
2154
2155
2156
2157
2158
2159
2159
2160
2161
2162
2163
2164
2165
2166
2167
2168
2169
2169
2170
2171
2172
2173
2174
2175
2176
2177
2178
2178
2179
2180
2181
2182
2183
2184
2185
2186
2187
2188
2188
2189
2190
2191
2192
2193
2194
2195
2196
2196
2197
2198
2199
2199
2200
2201
2202
2203
2204
2205
2206
2207
2208
2209
2209
2210
2211
2212
2213
2
```

```
62             leftScreen: {
63                 rejoinButtonEnabled: true,
64             },
65         };
66     );
67     const meeting = new VideoSDKMeeting();
68     meeting.init(config);
69 );
70 );
71 script.src =
72     "https://sdk.videosdk.live/rtc-js-prebuilt/0.3.31/rtc-js-
73     prebuilt.js";
74 document.getElementsByTagName("head")[0].appendChild(script);
75
76 </script>
77 </body>
78
79 </html>
80
```

Chat Bot

ChatBot.html

```
1  <!DOCTYPE html>
2  <html lang="en" dir="ltr">
3  <head>
4      <meta charset="utf-8">
5      <title>Chatbot in JavaScript | CodingNepal</title>
6      <link rel="stylesheet" href="ChatBot.css">
7      <meta name="viewport" content="width=device-width, initial-scale=1.0">
8      <link rel="stylesheet" href="https://fonts.googleapis.com/css2?family=Material+Symbols+Outlined:opsz,wght,FILL,GRAD@48,400,0,0" />
9      <link rel="stylesheet" href="https://fonts.googleapis.com/css2?family=Material+Symbols+Rounded:opsz,wght,FILL,GRAD@48,400,1,0" />
10     <script src="script.js" defer></script>
11 </head>
12 <body>
13     <h1>Remote Health Monitoring System Chat Bot</h1>
14     <button class="chatbot-toggler">
15         <span class="material-symbols-rounded">mode_comment</span>
16         <span class="material-symbols-outlined">close</span>
17     </button>
18     <div class="chatbot">
19         <header>
20             <h2>Chatbot</h2>
21             <span class="close-btn material-symbols-outlined">close</span>
22         </header>
23         <ul class="chatbox">
24             <li class="chat incoming">
25                 <span class="material-symbols-outlined">smart_toy</span>
26                 <p>Hi there !<br>How can I help you today?</p>
27             </li>
28         </ul>
29         <div class="chat-input">
30             <textarea placeholder="Enter a message..." spellcheck="false" required>
31             </textarea>
32             <span id="send-btn" class="material-symbols-rounded">send</span>
33         </div>
34     </div>
35 </body>
36 </html>
```

ChatBot.css

```
1 @import url('https://fonts.googleapis.com/css2?
family=Poppins:wght@400;500;600&display=swap');
2 ▼ *
3   margin: 0;
4   padding: 0;
5   box-sizing: border-box;
6   font-family: "Poppins", sans-serif;
7 }
8 ▼ body {
9   background: #E3F2FD;
10 }
11 ▼ .chatbot-toggler {
12   position: fixed;
13   bottom: 30px;
14   right: 35px;
15   outline: none;
16   border: none;
17   height: 50px;
18   width: 50px;
19   display: flex;
20   cursor: pointer;
21   align-items: center;
22   justify-content: center;
23   border-radius: 50%;
24   background: #724ae8;
25   transition: all 0.2s ease;
26 }
27 ▼ body.show-chatbot .chatbot-toggler {
28   transform: rotate(90deg);
29 }
30 ▼ .chatbot-toggler span {
31   color: #fff;
32   position: absolute;
33 }
34 .chatbot-toggler span:last-child,
35 ▼ body.show-chatbot .chatbot-toggler span:first-child {
36   opacity: 0;
37 }
38 ▼ body.show-chatbot .chatbot-toggler span:last-child {
39   opacity: 1;
40 }
41 ▼ .chatbot {
42   position: fixed;
43   right: 35px;
44   bottom: 90px;
45   width: 420px;
46   background: #fff;
47   border-radius: 15px;
48   overflow: hidden;
49   opacity: 0;
50   pointer-events: none;
51   transform: scale(0.5);
52   transform-origin: bottom right;
53   box-shadow: 0 0 128px 0 rgba(0,0,0,0.1),
54               0 32px 64px -48px rgba(0,0,0,0.5);
55   transition: all 0.1s ease;
56 }
57 ▼ body.show-chatbot .chatbot {
58   opacity: 1;
59   pointer-events: auto;
60   transform: scale(1);
61 }
```

```
62 ▼ .chatbot header {
63   padding: 16px 0;
64   position: relative;
65   text-align: center;
66   color: #fff;
67   background: #724ae8;
68   box-shadow: 0 2px 10px rgba(0,0,0,0.1);
69 }
70 ▼ .chatbot header span {
71   position: absolute;
72   right: 15px;
73   top: 50%;
74   display: none;
75   cursor: pointer;
76   transform: translateY(-50%);
77 }
78 ▼ header h2 {
79   font-size: 1.4rem;
80 }
81 ▼ .chatbot .chatbox {
82   overflow-y: auto;
83   height: 510px;
84   padding: 30px 20px 100px;
85 }
86 ▼ .chatbot :where(.chatbox, textarea)::webkit-scrollbar {
87   width: 6px;
88 }
89 ▼ .chatbot :where(.chatbox, textarea)::webkit-scrollbar-track {
90   background: #fff;
91   border-radius: 25px;
92 }
93 ▼ .chatbot :where(.chatbox, textarea)::webkit-scrollbar-thumb {
94   background: #ccc;
95   border-radius: 25px;
96 }
97 ▼ .chatbox .chat {
98   display: flex;
99   list-style: none;
100 }
101 ▼ .chatbox .outgoing {
102   margin: 20px 0;
103   justify-content: flex-end;
104 }
105 ▼ .chatbox .incoming span {
106   width: 32px;
107   height: 32px;
108   color: #fff;
109   cursor: default;
110   text-align: center;
111   line-height: 32px;
112   align-self: flex-end;
113   background: #724ae8;
114   border-radius: 4px;
115   margin: 0 10px 7px 0;
116 }
```

```
116  }
117 ▼ .chatbox .chat p {
118   white-space: pre-wrap;
119   padding: 12px 16px;
120   border-radius: 10px 10px 0 10px;
121   max-width: 75%;
122   color: #fff;
123   font-size: 0.95rem;
124   background: #724ae8;
125 }
126 ▼ .chatbox .incoming p {
127   border-radius: 10px 10px 10px 0;
128 }
129 ▼ .chatbox .chat p.error {
130   color: #721c24;
131   background: #f8d7da;
132 }
133 ▼ .chatbox .incoming p {
134   color: #000;
135   background: #f2f2f2;
136 }
137 ▼ .chatbot .chat-input {
138   display: flex;
139   gap: 5px;
140   position: absolute;
141   bottom: 0;
142   width: 100%;
143   background: #fff;
144   padding: 3px 20px;
145   border-top: 1px solid #ddd;
146 }
147 ▼ .chat-input textarea {
148   height: 55px;
149   width: 100%;
150   border: none;
151   outline: none;
152   resize: none;
153   max-height: 180px;
154   padding: 15px 15px 15px 0;
155   font-size: 0.95rem;
156 }
157 ▼ .chat-input span {
158   align-self: flex-end;
159   color: #724ae8;
160   cursor: pointer;
161   height: 55px;
162   display: flex;
163   align-items: center;
164   visibility: hidden;
165   font-size: 1.35rem;
166 }
167 ▼ .chat-input textarea:valid ~ span {
168   visibility: visible;
169 }
```

```
170
171 ▼ @media (max-width: 490px) {
172   .chatbot-toggler {
173     right: 20px;
174     bottom: 20px;
175   }
176   .chatbot {
177     right: 0;
178     bottom: 0;
179     height: 100%;
180     border-radius: 0;
181     width: 100%;
182   }
183   .chatbot .chatbox {
184     height: 90%;
185     padding: 25px 15px 100px;
186   }
187   .chatbot .chat-input {
188     padding: 5px 15px;
189   }
190   .chatbot header span {
191     display: block;
192   }
193 }
```

Script.js

```
1  const chatbotToggler = document.querySelector(".chatbot-toggler");
2  const closeBtn = document.querySelector(".close-btn");
3  const chatbox = document.querySelector(".chatbox");
4  const chatInput = document.querySelector(".chat-input textarea");
5  const sendChatBtn = document.querySelector(".chat-input span");
6
7  let userMessage = null;
8  const API_KEY = "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";
9  const inputInitHeight = chatInput.scrollHeight;
10
11 ▼ const createChatLi = (message, className) => {
12    const chatLi = document.createElement("li");
13    chatLi.classList.add("chat", `${className}`);
14    let chatContent = className === "outgoing" ? `<p></p>` : `<span
15      class="material-symbols-outlined">smart_toy</span><p></p>`;
16    chatLi.innerHTML = chatContent;
17    chatLi.querySelector("p").textContent = message;
18    return chatLi;
19  }
20
21 ▼ const generateResponse = (chatElement) => {
22   const API_URL = "https://api.openai.com/v1/chat/completions";
23   const messageElement = chatElement.querySelector("p");
24
25   const requestOptions = {
26     method: "POST",
27     headers: {
28       "Content-Type": "application/json",
29       "Authorization": `Bearer ${API_KEY}`
30     },
31     body: JSON.stringify({
32       model: "gpt-3.5-turbo",
33       messages: [{role: "user", content: userMessage}],
34     })
35
36   fetch(API_URL, requestOptions).then(res => res.json()).then(data => {
37     messageElement.textContent = data.choices[0].message.content.trim();
38   }).catch(() => {
39     messageElement.classList.add("error");
40     messageElement.textContent = "Oops! Something went wrong. Please try
41     again.";
42   }).finally(() => chatbox.scrollTo(0, chatbox.scrollHeight));
43 }
44
45 ▼ const handleChat = () => {
46   userMessage = chatInput.value.trim();
47   if(!userMessage) return;
48
49   chatInput.value = "";
50   chatInput.style.height = `${inputInitHeight}px`;
51
52   chatbox.appendChild(createChatLi(userMessage, "outgoing"));
53   chatbox.scrollTo(0, chatbox.scrollHeight);
54
55   setTimeout(() => {
56     const incomingChatLi = createChatLi("Thinking...", "incoming");
57     chatbox.appendChild(incomingChatLi);
58     chatbox.scrollTo(0, chatbox.scrollHeight);
59   }, 600);
60 }
```

```
61
62 ˇ chatInput.addEventListener("input", () => {
63      chatInput.style.height = `${inputInitHeight}px`;
64      chatInput.style.height = `${chatInput.scrollHeight}px`;
65  });
66
67 ˇ chatInput.addEventListener("keydown", (e) => {
68 ˇ     if(e.key === "Enter" && !e.shiftKey && window.innerWidth > 800) {
69         e.preventDefault();
70         handleChat();
71     }
72 });
73
74 sendChatBtn.addEventListener("click", handleChat);
75 closeBtn.addEventListener("click", () => document.body.classList.remove("show-
chatbot"));
76 chatbotToggler.addEventListener("click", () =>
document.body.classList.toggle("show-chatbot"));
```

18. User guide and installation guide

Step 1: Download and zip the file




Step 2: Put the zipped file into the path


X:\xampp\htdocs


名称	修改日期	类型	大小
dashboard	2023/4/6 0:34	文件夹	
FYP web	2024/4/19 22:36	文件夹	
test1	2024/4/13 23:12	文件夹	
xampp	2023/4/6 0:34	文件夹	
applications	2022/6/16 0:07	Firefox HTML D...	4 KB
bitnami	2022/6/16 0:07	CSS Document	1 KB
favicon	2015/7/16 23:32	ICO 文件	31 KB
index	2015/7/16 23:32	PHP 文件	1 KB

Step 3: Open the Xampp and Turn on the Apache and MySQL




Step4: Open the Admin and insert the mysql file






Step 5: Go to the browser and paste the link: <http://127.0.0.1/FYP web/>






Step6: Click the button 'Sign Up New' to sign up



Choose your character. If you are a public user follow step 7, otherwise follow step 12.

Step7: Login or sign up a new account (For Public)





Step8: Make Appointment (For Public)

The screenshot shows a web-based application interface. At the top, there is a navigation bar with links: Home, Appointment, HealthyTest, Meeting, Chat, Healthy Information, Profile, and Logout. Below the navigation bar is a green header bar labeled "Doctor List". A large blue arrow points downwards from the top of the page towards the "Doctor List" table. The table has columns: Id, Name, Types, Tel, Email, and Appointment. There are four rows of data:

Id	Name	Types	Tel	Email	Appointment
1	Mike Lee	General medical practitioner	51212456	Mike12@gmail.com	<input type="button" value=""/>
2	William Hunt	General medical practitioner	61214456	William12@gmail.com	<input type="button" value=""/>
3	Denise Dunn	Cardiologist	62215477	Denise22@gmail.com	<input type="button" value=""/>
4	som	General medical practitioner	51712126	s123@gmail.com	<input type="button" value=""/>



1. Choose the doctor.
2. Click the related row appointment button to open the form


The screenshot shows the same application interface as the previous one, but with a modal dialog box overlaid on the "Doctor List" table. The modal is titled "Appointment" and contains the following fields:

- Name: Isla Jones
- Date: A date input field with a calendar icon, currently showing "yyyy / mm / dd".
- Time: A time input field showing "-- : --".
- Doctor Name: Mike Lee
- Visit Method: A dropdown menu labeled "Select method".
- Submit button (yellow)
- Clear button (pink)



3. Fill the form
4. Click submit button

Step9: Click the button chat to communicate with doctor (Public)






Step10: Click the button meeting to communicate with doctor (Public)



Step11: Click the Healthy Information to Search the information and use the chat bot



A screenshot of the 'Healthy Info' page. A red arrow points to the 'Neurological Diseases' section in the sidebar. The main content area shows a search bar with 'Enter disease...' and a 'search' button. A result for 'Polio' is displayed, with a red arrow pointing to the title 'Polio'. The sidebar also lists other categories and diseases: Parkinson's Disease, Endocrine Diseases, Diabetes, Cardiovascular Diseases, Heart Disease, Mental and Psychological Diseases, Depression, Mental Illness, Anorexia, Respiratory Diseases, and Asthma.

Chatbot

Hi there 🤖
How can I help you
today?



改善血滿足身體的需
要，可能導致心臟病



是心臟病的主要風


Enter a message...

下，可能需要進行血



Step12: Login or sign up a new account (For Doctor)





Step13: Management appointment (Doctor)

Step13: management appointment (doctor)

A screenshot of a web-based appointment management system. At the top, there's a navigation bar with icons for Home, Appointment Record, Chat, Meeting, Sensor Record, and Profile. Below the navigation is a table titled "Appointment List". The table has columns for Id, Name, Date, Time, Visit Method, Doctor Name, Edit, and Delete. A single row is shown with the following data: Id: 4, Name: Isla Jones, Date: 2024-04-22, Time: 04:08:00.000000, Visit Method: online, Doctor Name: Mike Lee. To the right of the row are green "Edit" and red "Delete" buttons. Two arrows point downwards from the text "Doctor can edit the related row record or delete" towards these buttons.

Appointment List							
Id	Name	Date	Time	Visit Method	Doctor Name	Edit	Delete
4	Isla Jones	2024-04-22	04:08:00.000000	online	Mike Lee		

Doctor can edit the related row record or delete

A screenshot of a "Edit Appointment" modal dialog. At the top, it says "Edit Appointment" with a close button. The dialog contains the following fields: ID: 4, Name: Isla Jones, Date: 2024/04/22, Time: 04:08, Visit Method: online, Doctor Name: Mike Lee. At the bottom are two buttons: a yellow "Submit" button and a red "Clear" button.

Doctor Name	Edit	Delete
Mike Lee		

Step14: Sensor Healthy test (Public)

Sensor data display

Sensor Data Dashboard

Temperature Sensor
Temperature: 36.7 °C

Heart Rate Sensor
Heart Rate: 61.68333333333333 BPM

Human Resistance Sensor
Resistance: 610 Ohms

History Record				
ID	Date & Time	Temperature	HeartRate	GSR
1	2024-04-06 17:00:25	36.7	80	610
2	2024-04-07 17:00:33	36.7	80	610
3	2024-04-07 17:14:20	36.7	80	610
4	2024-04-07 17:16:55	36.7	80	610
5	2024-04-07 17:20:56	36.7	80	610
6	2024-04-07 17:21:59	36.7	80	610

Upload to history record table.

Sensor Data Dashboard

Temperature Sensor
Temperature: 36.7 °C

Heart Rate Sensor
Heart Rate: 61.68333333333333 BPM

Human Resistance Sensor
Resistance: 610 Ohms

History Record				
ID	Date & Time	Temperature	HeartRate	GSR
1	2024-04-06 17:00:25	36.7	80	610
2	2024-04-07 17:00:33	36.7	80	610
3	2024-04-07 17:14:20	36.7	80	610
4	2024-04-07 17:16:55	36.7	80	610
5	2024-04-07 17:20:56	36.7	80	610
6	2024-04-07 17:21:59	36.7	80	610

Step15: View Sensor Record (Doctor)

Remote health monitoring system

- [Home](#)
- [Appointment Record](#)
- [Chat](#)
- [Meeting](#)
- [Sensor Record](#)
- [Profile](#)

Patient List

Sensor Record

ID	Name	Email	Tel
1	Isla Jones	Isla2j@gmail.com	55678899
2	Juan Owens	Juow11@gmail.com	65446322
3	TAM	s123@gmail.com	51712126
4	ssaa	s123@gmail.com	51234561

Remote health monitoring system

Patient List

Sensor Record

Search by patient name...

ID	Name	Date & Time	Temperature	Heart Rate	GSR
1	Isla Jones	2024-04-06 17:00:25	36.7	80	610
2	Isla Jones	2024-04-07 17:00:33	36.7	80	610
3	Isla Jones	2024-04-07 17:14:20	36.7	80	610
4	Isla Jones	2024-04-07 17:16:55	36.7	80	610
5	Isla Jones	2024-04-07 17:20:56	36.7	80	610
6	Isla Jones	2024-04-07 17:21:59	36.7	80	610
7	Isla Jones	2024-04-07 18:46:33	36.7	80	610
8	Juan Owens	2024-04-07 21:08:52	36.9	76	608

Remote health monitoring system

Patient List

Sensor Record

ju

ID	Name	Date & Time	Temperature	Heart Rate	GSR
8	Juan Owens	2024-04-07 21:08:52	36.9	76	608